**Lipidic Control
IVD
Store at 2 - 8°C**
PRODUCT CHARACTERISTICS

LIPIDIC CONTROL is a human lyophilised control serum. Is intended to monitor the accuracy and precision of Cholesterol, HDL cholesterol, LDL cholesterol, Phospholipids and Triglycerides.

REAGENTS

Human serum. Biological additives. Bacteriostatic agents. The concentration of the components are lot-specific. See chart below for the values of ingredients.

PRECAUTIONS

Components from human origin have been tested and found negative for the presence of HBsAg, HCV and antibody to HIV (1/2). However handle cautiously as potentially infectious.

CALIBRATION

Information available on request.

PREPARATION

- Open the bottle carefully, to avoid loss of lyophilized.
- Reconstitute the lyophilised vial with 1 mL distilled water.
- Bring to room temperature and mix thoroughly, avoiding foam forming. Wait for 30 minutes before use.

Improper handling and/or storage can affect results. Inaccurate reconstitution and errors in assay technique can cause erroneous results.

STORAGE AND STABILITY

The Control serum is stable until the expiration date indicated on the label when stored tightly closed at 2-8°C and contaminations are prevented during its use. Do not use reagents over the expiration date or if there is visible evidence of microbial growth.

Store tightly capped when not use.

- After reconstitution is stable for:

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