D-DIMER CONTROL

### D-Dimer / Dímero-D Control Turbidimetry / Turbidimetría

D-DIMER CONTROL

### D-Dimer / Dímero-D Control Turbidimetry / Turbidimetría

#### Quantitative determination of D-Dimer

IVD  
Store at 2 - 8°C

#### PRODUCT CHARACTERISTICS

The control is a lyophilised human serum used to evaluate the precision and accuracy of D-Dimer determinations in human plasma or serum.

#### REAGENTS

DIMER-D CTRL	Human serum. Contains ProClin 300, as a preservative. The assigned values of D-Dimer are shown on the below table
--------------	--

#### PRECAUTIONS

H317: Contains: 5-chloro-2-methyl-4-isothiazolin-3-one and 2-methyl-2H-isothiazol-3-one (3:1), reaction mass of: 5-chloro-2-methyl-4-isothiazolin-3-one and 2-methyl-4-isothiazolin-3-one (3:1) (Proclin 300). May produce an allergic reaction. Follow the precautionary statements given in MSDS and label of the product.

Components from human origin have been tested and found negative for the presence of HBsAg, HCV and antibody to HIV (1/2). However handle cautiously as potentially infectious.

#### PREPARATION

Reconstitute with 1.0 mL of distilled water. Mix thoroughly, avoiding foam forming. Bring at room temperature for about 10 minutes before use.

#### STORAGE AND STABILITY

The control is stable until the expiration date on the label when stored tightly closed at 2-8°C and contaminations are prevented during their use. Do not use reagents over the expiration date.

After reconstitution, stable for 14 days at 2-8°C or 1 month at less than -80°C. Freeze and thaw of the calibrator solution can be permitted only one time for the measurement.

#### PROCEDURE

This control should be used in the D-Dimer reagent.

PRODUCT/PRODUCTO	VALUE/VALOR	RANGE/RANGO
D-Dimer LOW	3.3 ug/mL	1.8 – 4.8
D-Dimer HIGH	10.1 ug/mL	8.1– 12.1

Actual results depend on many factors, including lot number, type of reagent, and instrument. Ranges must be determined in each laboratory with each change of lot number, reagent, or instrument.

This Control value sheet is applicable to sublots. Sequential alphabetical letter (e.g. A, B, C etc.) following the lot number.

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