
**Quantitative determination of alanine aminotransferase
GPT (ALT)
IVD**

Store at 2-8°C

PRINCIPLE OF THE METHOD

Alanine aminotransferase (ALT) or Glutamate pyruvate transaminase (GPT) catalyses the reversible transfer of an amino group from alanine to α -ketoglutarate forming glutamate and piruvate. The piruvate produced is reduced to lactate by lactate dehydrogenase (LDH) and NADH:

The rate of decrease in concentration of NADH, measured photometrically, is proportional to the catalytic concentration of ALT present in the sample¹.

CLINICAL SIGNIFICANCE

The ALT is a cellular enzyme, found in highest concentration in liver and kidney.

High levels are observed in hepatic disease like hepatitis, diseases of muscles and traumatism; its better application is in the diagnosis of the diseases of the liver.

When they are used in conjunction with AST aid in the diagnosis of infarcts in the myocardium, since the value of the ALT stays within the normal limits in the presence of elevated levels of AST^{1,4,5}.

Clinical diagnosis should not be made on a single test result; it should integrate clinical and other laboratory data.

REAGENTS

R 1	TRIS pH 7.8	100 mmol/L
Buffer	Lactate dehydrogenase (LDH)	1200 U/L
	L-Alanine	500 mmol/L
R 2	NADH	0,18 mmol/L
Substrate	α -Ketoglutarate	15 mmol/L

PRECAUTIONS

R1: H290-May be corrosive to metals.

Follow the precautionary statements given in MSDS and label of the product.

PREPARATION

All the reagents are ready to use.

STORAGE AND STABILITY

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C, protected from light and contaminations prevented during their use.

Do not use reagents over the expiration date.

Signs of reagent deterioration:

- Presence of particles and turbidity.
- Blank absorbance (A) at 340 nm < 1.00.

ADDITIONAL EQUIPMENT

- SPIN640 / SPIN640Plus Autoanalyzer.
- General laboratory equipment.

SAMPLES

Serum or plasma¹: Stability 7 days at 2-8°C.

QUALITY CONTROL

Control sera are recommended to monitor the performance of assay procedures: SPINROL H Normal and Pathologic (Ref. 1002120 and 1002210).

If control values are found outside the defined range, check the instrument, reagents and technique for problems.

Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

BARCODED REAGENTS LOAD MUST BE PRECEDED OF A SPINREACT "DATABASE" COPY INTO THE ANALYZER SOFTWARE. IT IS AVAILABLE UNDER REQUEST TO SPINREACT.

SPIN640 APPLICATION

TEST INFORMATION		REAGENT VOLUME	
Nº	**	Vol. R1	240
Test	GPT	Vol. R2	60
Full Name	GPT	Vol. R3	
Standard nº		Vol. R4	
SAMPLE VOLUME		RESULT SETUP	
Vol. Sample Stand.	30	Decimal	1 Slope 1
Vol. Sample Increas.		Unit	U/L Inter. 0
Vol. Sample Dec			
REACTION PARAMETERS			
Reac. Type	Kinetic	Direction	Decrease
Pri. Wave.	340	Reagent Blank	0-0
Sec. Wave.		React. Time	50-70

SPIN640Plus APPLICATION

EDIT PARAMETERS			
Test	GPT	No.	**
Full name	GPT	Print name	GPT
Reac. Type	Kinetic	Direction	Derease
Pri. Wave.	340	Sec. Wave.	
Unit	U/L	Decimal	0.1
Reagent Blank	0 - 0	React. Time	57 - 77
Vol. Sample	30 ul	R1	240 ul
Increased		R2	60 ul
Decreased		R3	
Sample blank		R4	

The Calibration is stable until **35 days**. After this period the Calibration must be performed again in order to obtain good results.

REFERENCE VALUES^{4,5}

- Men 40 U/L
- Women 32 U/L

Normal newborns have been reported to show a reference range of up to double the adult, attributed to the neonate's hepatocytes. These values decline to adult levels by approximately 3 months of age.

These values are for orientation purpose; each laboratory should establish its own reference range.

PERFORMANCE CHARACTERISTICS

Measuring range: From *detection limit* of 0 U/L to *linearity limit* of 400 U/L. If the results obtained were greater than linearity limit, dilute the sample 1/10 with NaCl 9 g/L and multiply the result by 10.

Precision:

	Intra-assay (n=20)		Inter-assay (n=20)	
Mean (U/L)	42,0	116	41,1	115
SD	0,47	0,42	0,76	1,61
CV (%)	1,11	0,36	1,85	1,40

Sensitivity: 1 U/L = 0,00052 ΔA / min.

Accuracy: Results obtained using SPINREACT reagents (y) did not show systematic differences when compared with other commercial reagents (x).

The results obtained using 50 samples were the following:

Correlation coefficient (r)²: 0,99597.

Regression equation: y=1,1209x + 1,390.

The results of the performance characteristics depend on the analyzer used.

BIBLIOGRAPHY

1. Murray R. Alanine aminotransferase. Kaplan A et al. Clin Chem The C.V. Mosby Co. St Louis. Toronto. Princeton 1984; 1088-1090.
2. Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.
3. Young DS. Effects of disease on Clinical Lab. Tests, 4th ed AACC 2001.
4. Burtis A et al. Tietz Textbook of Clinical Chemistry, 3rd ed AACC 1999.
5. Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed AACC 1995.

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