
**Quantitative determination of lactate dehydrogenase (LDH)**
**IVD**

Store at 2-8°C

**PRINCIPLE OF THE METHOD**

Lactate dehydrogenase (LDH) catalyses the reduction of pyruvate by NADH, according the following reaction:

 The rate of decrease in concentration of NADPH, measured photometrically, is proportional to the catalytic concentration of LDH present in the sample<sup>1</sup>.

**CLINICAL SIGNIFICANCE**

Lactate dehydrogenase (LDH) is an enzyme with wide tissue distribution in the body.

The higher concentrations of LDH are found in liver, heart, kidney, skeletal muscle and erythrocytes.

 Increased levels of the enzyme are found in serum in liver disease, myocardial infarction, renal disease, muscular dystrophy and anemia<sup>1,4,5</sup>.

Clinical diagnosis should not be made on a single test result; it should integrate clinical and other laboratory data.

**REAGENTS**

<b>Reagent 1</b>	Imidazol	65 mmol/L
Buffer	Pyruvate	0,6 mmol/L
<b>Reagent 2</b>	NADH	0,18 mmol/L
Substrate		

**PREPARATION**

All the reagents are ready to use.

**STORAGE AND STABILITY**

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C, protected from light and contaminations prevented during their use.

Do not use reagents over the expiration date.

**Signs of reagent deterioration:**

- Presence of particles and turbidity.
- Blank absorbance (A) at 340 nm <1,00.

**ADDITIONAL EQUIPMENT**

- MINDRAY BS-120 / BS-200E Autoanalyzer.
- General laboratory equipment.

**SAMPLES**

 Serum<sup>1</sup>. Separated from cells as rapidly as possible. Do not use oxalates as anticoagulants since they inhibit the enzyme.

Do not use haemolysed samples. Stability: 2 days at 2-8°C.

**REFERENCE VALUES<sup>1</sup>**

25°C	30°C	37°C
120-240 U/L	160-320 U/L	230-460 U/L

These values are for orientation purpose; each laboratory should establish its own reference range.

**QUALITY CONTROL**

Control sera are recommended to monitor the performance of assay procedures: SPINTROL H Normal and Pathologic (Ref. 1002120 and 1002210).

If control values are found outside the defined range, check the instrument, reagents and technique for problems.

Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

**MINDRAY BS-120 / BS-200E APPLICATION**
**PARAMETERS**

Test	LDH / LDH	R1	240 / 240
Nº	**	R2	60 / 60
Full Name	LDH / LDH	Sample volume	5 / 5
Standard Nº		R1 Blank	
Reac. Type	Kinet / Kinet	Mixed Rgt Blank	
Pri. Wavelength	340 / 340	Linearity Range	5 U/L 1450 U/L
Sec. Wavelength		Linearity Limit	*
Direction	Decrea / Decrea	Substrate Limit	*
Reac. Time	4_14 / 4_14	Factor	*
Incuba. Time		Prozone check	*
Units	U/L / U/L	q1	q2
Precision	Interger / Interger	q3	q4
		PC	Abs

**CALIBRATION**

Rule	One-point Linear / Two-point Linear
Sensitivity	1 / 1
Replicates	2 / 2
Interval (days)	0 / 0
Difference Limit	
SD	
Blank Response	
Error Limit	
Correlation Coefficient	

*Blank parameter must be performed in order to get good results in CALIB screen from main menu. The blank calibration is stable until 35 days. After this period the blank parameter must be performed again in order to validate the calibration.*

**PERFORMANCE CHARACTERISTICS**
**Measuring range:** From detection limit of 3,42 U/L to linearity limit of 1600 U/L.

If the results obtained were greater than linearity limit, dilute the sample 1/10 with NaCl 9 g/L and multiply the result by 10.

**Precision:**

	Intra-assay (n=20)		Inter-assay (n=20)	
Mean (U/L)	400	785	392	773
SD	3,15	10,97	6,23	9,93
CV (%)	0,79	1,40	1,59	1,28

**Sensitivity:** 1 U/L = 0,00009 ΔA/min.

**Accuracy:** Results obtained using SPINREACT reagents (y) did not show systematic differences when compared with other commercial reagents (x).

The results obtained using 50 samples were the following:

 Correlation coefficient (r)<sup>2</sup>: 0,98382.

Regression equation: y= 0,8988x + 2,583.

The results of the performance characteristics depend on the analyzer used.

**BIBLIOGRAPHY**

1. Pesce A. Lactate dehydrogenase. Kaplan A et al. Clin Chem The C.V. Mosby Co. St Louis. Toronto. Princeton 1984; 1124-117, 438.
2. Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.
3. Young DS. Effects of disease on Clinical Lab. Tests, 4th ed AACC 2001.
4. Burtis A et al. Tietz Textbook of Clinical Chemistry, 3rd ed AACC 1999.
5. Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed AACC 1995.

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