
**Quantitative determination of alanine aminotransferase GPT (ALT) IVD**

Store at 2-8°C

**PRINCIPLE OF THE METHOD**

Alanine aminotransferase (ALT) or Glutamate pyruvate transaminase (GPT) catalyses the reversible transfer of an amino group from alanine to  $\alpha$ -ketoglutarate forming glutamate and pyruvate.

The pyruvate produced is reduced to lactate by lactate dehydrogenase (LDH) and NADH:

The rate of decrease in concentration of NADH, measured photometrically, is proportional to the catalytic concentration of ALT present in the sample<sup>1</sup>.

**CLINICAL SIGNIFICANCE**

The ALT is a cellular enzyme, found in highest concentration in liver and kidney.

High levels are observed in hepatic disease like hepatitis, diseases of muscles and traumatism; its better application is in the diagnosis of the diseases of the liver.

When they are used in conjunction with AST aid in the diagnosis of infarcts in the myocardium, since the value of the ALT stays within the normal limits in the presence of elevated levels of AST<sup>1,4,5</sup>.

Clinical diagnosis should not be made on a single test result; it should integrate clinical and other laboratory data.

**REAGENTS**

R 1 Buffer	TRIS pH 7,8	100 mmol/L
	Lactate dehydrogenase (LDH)	1200 U/L
R 2 Substrate	L-Alanine	500 mmol/L
	NADH	0,18 mmol/L
	$\alpha$ -Ketoglutarate	15 mmol/L

**PRECAUTIONS**

R1: H290-May be corrosive to metals.

Follow the precautionary statements given in MSDS and label of the product.

**PREPARATION**

All the reagents are ready to use.

**STORAGE AND STABILITY**

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C, protected from light and contaminations prevented during their use.

Do not use reagents over the expiration date.

**Signs of reagent deterioration:**

- Presence of particles and turbidity.
- Blank absorbance (A) at 340 nm < 1,00.

**ADDITIONAL EQUIPMENT**

- MINDRAY BS-120 / BS-200E Autoanalyzer.
- General laboratory equipment.

**SAMPLES**

Serum or plasma<sup>1</sup>: Stability 7 days at 2-8°C.

**Temperature conversion factors**

To correct results to other temperatures multiply by:

Assay temperature	Conversion factor to		
	25°C	30°C	37°C
25°C	1,00	1,32	1,82
30°C	0,76	1,00	1,39
37°C	0,55	0,72	1,00

**REFERENCE VALUES<sup>4,5</sup>**

	25°C	30°C	37°C
Men	up to 22 U/L	29 U/L	40 U/L
Women	up to 18 U/L	22 U/L	32 U/L

Normal newborns have been reported to show a reference range of up to double the adult, attributed to the neonate's hepatocytes. These values decline to adult levels by approximately 3 months of age.

These values are for orientation purpose; each laboratory should establish its own reference range.

**MINDRAY BS-120 / BS-200E APPLICATION**

PARAMETERS			
Test	GPT / GPT	R1	240 / 240
Nº	**	R2	60 / 60
Full Name	GPT / GPT	Sample volume	30 / 30
Standard Nº		R1 Blank	
Reac. Type	Kinet / Kinet	Mixed Rgt Blank	
Pri. Wavelength	340 / 340	Linearity Range	5 U/L 500 U/L
Sec. Wavelength		Linearity Limit	*
Direction	Decrea / Decrea	Substrate Limit	*
Reac. Time	3_10 / 3_10	Factor	*
Incuba. Time		Prozone check	*
Units	U/L / U/L	q1	q2
Precision	Interger / Interger	q3	q4
		PC	Abs
CALIBRATION			
Rule	One-point Linear / Two-point Linear		
Sensitivity	1 / 1		
Replicates	2 / 2		
Interval (days)	0 / 0		
Difference Limit			
SD			
Blank Response			
Error Limit			
Correlation Coefficient			

Blank parameter must be performed in order to get good results in CALIB screen from main menu. The blank calibration is stable until **35 days**. After this period the blank parameter must be performed again in order to validate the calibration.

**QUALITY CONTROL**

Control sera are recommended to monitor the performance of assay procedures: SPINTROL H Normal and Pathologic (Ref. 1002120 and 1002210).

If control values are found outside the defined range, check the instrument, reagents and technique for problems.

Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

**PERFORMANCE CHARACTERISTICS**

**Measuring range:** From detection limit of 0 U/L to linearity limit of 400 U/L.

If the results obtained were greater than linearity limit, dilute the sample 1/10 with NaCl 9 g/L and multiply the result by 10.

**Precision:**

	Intra-assay (n=20)		Inter-assay (n=20)	
Mean (U/L)	42,0	116	41,1	115
SD	0,47	0,42	0,76	1,61
CV (%)	1,11	0,36	1,85	1,40

**Sensitivity:** 1 U/L = 0,00052  $\Delta$ A / min.

**Accuracy:** Results obtained using SPINREACT reagents (y) did not show systematic differences when compared with other commercial reagents (x).

The results obtained using 50 samples were the following:

Correlation coefficient (r)<sup>2</sup>: 0,99597.

Regression equation: y=1,1209x + 1,390.

The results of the performance characteristics depend on the analyzer used.

**BIBLIOGRAPHY**

- Murray R. Alanine aminotransferase. Kaplan A et al. Clin Chem The C.V. Mosby Co. St Louis. Toronto. Princeton 1984; 1088-1090.
- Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.
- Young DS. Effects of disease on Clinical Lab. Tests, 4th ed AACC 2001.
- Burtis A et al. Tietz Textbook of Clinical Chemistry, 3rd ed AACC 1999.
- Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed AACC 1995.

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