
Quantitative determination of gamma-glutamyl transferase (γ-GT)
IVD

Store at 2-8°C

PRINCIPLE OF THE METHOD

Gamma-glutamyl transferase (γ-GT) catalyses the transfer of γ-glutamyl group from γ-glutamyl-p-nitroanilide to acceptor glycylglycine, according to the following reaction:

The rate of 2-nitro-5-aminobenzoic acid formation, measured photometrically, is proportional to the catalytic concentration of γ-GT present in the sample^{1,2}.

CLINICAL SIGNIFICANCE

Gamma-glutamyl transferase (γ-GT) is a cellular enzyme with wide tissue distribution in the body, primarily in the kidney, pancreas, liver and prostate.

Measurements of gamma-glutamyl transferase (γ-GT) activity are used in the diagnosis and treatment of hepatobiliary diseases such as biliary obstruction, cirrhosis or liver tumours^{1,2,5,6}.

Clinical diagnosis should not be made on a single test result; it should integrate clinical and other laboratory data.

REAGENTS

R 1	TRIS pH 8.6	100 mmol/L
Buffer	Glycylglycine	100 mmol/L
R 2	L-γ-glutamyl-3-carboxy-4-nitroanilide	3 mmol/L
Substrate		

PREPARATION

All the reagents are ready to use.

STORAGE AND STABILITY

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C, protected from light and contaminations prevented during their use.

Do not use reagents over the expiration date.

Signs of reagent deterioration:

- Presence of particles and turbidity.
- Blank absorbance (A) at 405 nm ≥ 1,20.

ADDITIONAL EQUIPMENT

- MINDRAY BS-120 / BS-200E Autoanalyzer.
- General laboratory equipment.

SAMPLES

Serum¹. γ-GT is stable for at least 3 days at 2-8°C, 8 hours at 15-25°C and 1 month at -20°C.

REFERENCE VALUES¹

	25°C	30°C	37°C
Women	4-18 U/L	5-25 U/L	7-32 U/L
Men	6-28 U/L	8-38 U/L	11-50 U/L

These values are for orientation purpose; each laboratory should establish its own reference range.

QUALITY CONTROL

Control sera are recommended to monitor the performance of assay procedures: SPINTROL H Normal and Pathologic (Ref. 1002120 and 1002210).

If control values are found outside the defined range, check the instrument, reagents and technique for problems.

Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

MINDRAY BS-120 / BS-200E APPLICATION

PARAMETERS			
Test	GGT / GGT	R1	240 / 240
Nº	**	R2	60 / 60
Full Name	GGT / GGT	Sample volume	30 / 30
Standard Nº		R1 Blank	
Reac. Type	Kinet / Kinet	Mixed Rgt Blank	
Pri. Wavelength	405 / 412	Linearity Range	2 U/L 250 U/L
Sec. Wavelength		Linearity Limit	*
Direction	Increase / Increase	Substrate Limit	*
Reac. Time	3_10 / 3_10	Factor	*
Incuba. Time		Prozone check	*
Units	U/L / U/L	q1	q2
Precision	Interger / Interger	q3	q4
		PC	Abs
CALIBRATION			
Rule	One-point Linear / Two-point Linear		
Sensitivity	1 / 1		
Replicates	2 / 2		
Interval (days)	0 / 0		
Difference Limit			
SD			
Blank Response			
Error Limit			
Correlation Coefficient			

Blank parameter must be performed in order to get good results in CALIB screen from main menu. The blank calibration is stable until **35 days**. After this period the blank parameter must be performed again in order to validate the calibration.

PERFORMANCE CHARACTERISTICS

Measuring range: From detection limit of 2 U/L to linearity limit of 300 U/L.

If the results obtained were greater than linearity limit, dilute the sample 1/10 with NaCl 9 g/L and multiply the result by 10.

Precision:

	Intra-assay (n=20)		Inter-assay (n=20)	
	Mean (U/L)	SD	Mean (U/L)	SD
Mean (U/L)	38,3	190	40,1	198
SD	0,39	0,53	0,82	2,30
CV (%)	1,03	0,28	2,05	1,16

Sensitivity: 1 U/L = 0,0008 ΔA/min.

Accuracy: Results obtained using SPINREACT reagents (y) did not show systematic differences when compared with other commercial reagents (x).

The results obtained using 50 samples were the following:

Correlation coefficient (r)²: 0,99990.

Regression equation: y= 1,334x - 1,493.

The results of the performance characteristics depend on the analyzer used.

BIBLIOGRAPHY

1. Gendler S. γ-GT. Kaplan A et al. Clin Chem The C.V. Mosby Co. St Louis. Toronto. Princeton 1984; 1120-1123.
2. Persijn J P et al. J Clin Chem Clin Biochem 1976; (14) 9: 421-427.
3. Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.
4. Young DS. Effects of disease on Clinical Lab. Tests, 4th ed AACC 2001.
5. Burtis A et al. Tietz Textbook of Clinical Chemistry, 3rd ed AACC 1999.
6. Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed AACC 1995.

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