
**Quantitative determination of chloride ion IVD**

Store at 2-8°C

**PRINCIPLE OF THE METHOD**

The quantitative displacement of thiocyanate by chloride from mercuric thiocyanate and subsequent formation of a red ferric thiocyanate complex is measured colorimetrically:

The intensity of the color formed is proportional to the chloride ion concentration in the sample<sup>1,2,3,4</sup>.

**CLINICAL SIGNIFICANCE**

It is important clinically the determination of chloride due regulation of osmotic pressure of extra cellular fluid and to its significant role in acid-base balance. Increases in chloride ion concentration may be found in severe dehydration, excessive intake of chloride, severe renal tubular damage and in patients with cystic fibrosis.

Decrease in chloride ion concentration may be found in metabolic acidosis, loss from prolonged vomiting and chronic pyelonephritis<sup>2,7,8</sup>. Clinical diagnosis should not be made on a single test result; it should integrate clinical and other laboratory data.

**REAGENTS**

<b>R</b> Thiocyanate-Hg	Mercuric thiocyanate	4 mmol/L
	Ferric nitrate	40 mmol/L
	Mercuric nitrate	2 mmol/L
	Nitric acid	45 mmol/L

**PRECAUTIONS**

R: H314-Causes severe skin burns and eye damage.  
Follow the precautionary statements given in MSDS and label of the product.

**PREPARATION**

The reagent is ready to use.

**STORAGE AND STABILITY**

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C, protected from light and contaminations prevented during their use. Do not use reagents over the expiration date.

**Signs of reagent deterioration:**

- Presence of particles and turbidity.
- Blank absorbance (A) at 505 nm  $\geq 0,15$ .

**ADDITIONAL EQUIPMENT**

- MINDRAY BS-120 / BS-200E Autoanalyzer.
- General laboratory equipment (Note 1, 2, 3).

**SAMPLES**

- Serum or plasma free of haemolysis and separated from cells as rapidly as possible. Anticoagulants such as oxalate or EDTA are not acceptable they will interfere with results.
- Urine<sup>1</sup>: Collect 24-hour urine specimen in chloride free containers. Dilute a sample 1/2 in distilled water. Mix. Multiply results by 2 (dilution factor).

Stability of the sample: Ion chloride is stable 1 week at room temperature (15-25°C), in refrigerator (2-8°C) or frozen (-20°C) temperatures.

**REFERENCE VALUES<sup>1</sup>**

Serum or plasma:	95 - 115 mmol/L
Urine:	110 - 250 mmol/24h

These values are for orientation purpose; each laboratory should establish its own reference range.

**QUALITY CONTROL**

Control sera and calibrators are recommended to monitor the performance of assay procedures: SPINTROL H Calibrator, SPINTROL H Normal and Pathologic (Ref. 1002011, 1002120 and 1002210).

If control values are found outside the defined range, check the instrument, reagents and technique for problems.

Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

**MINDRAY BS-120 / BS-200E APPLICATION**

<u>PARAMETERS</u>			
Test	CL / CL	R1	300 / 300
Nº	**	R2	*
Full Name	CL / CL	Sample volume	3 / 3
Standard Nº		R1 Blank	
Reac. Type	Endpoint / Endpoint	Mixed Rgt Blank	
Pri. Wavelength	510 / 505	Linearity Range	0 mmol/L 130 mmol/L
Sec. Wavelength		Linearity Limit	*
Direction	Increase / Increase	Substrate Limit	*
Reac. Time	0_17 / 0_17	Factor	*
Incuba. Time		Prozone check	*
Units	mmol/L / mmol/L	q1	q2
Precision	Interger / Interger	q3	q4
		PC	Abs
<u>CALIBRATION (Cal + Rgt Blk)</u>			
Rule	One-point Linear / Two-point Linear		
Sensitivity	1 / 1		
Replicates	2 / 2		
Interval (days)	0 / 0		
Difference Limit			
SD			
Blank Response			
Error Limit			
Correlation Coefficient			

Blank parameter must be performed in order to get good results in CALIB screen from main menu. The blank calibration is stable until **35 days**. After this period the blank parameter must be performed again in order to validate the calibration.

**PERFORMANCE CHARACTERISTICS**

**Measuring range:** From *detection limit* of 0,454 mmol/L to *linearity limit* of 190 mmol/L.

If the results obtained were greater than linearity limit, dilute the sample 1/2 with distilled water and multiply the result by 2.

**Precision:**

Mean (mmol/L)	Intra-assay (n=20)		Inter-assay (n=20)	
	SD	CV (%)	SD	CV (%)
84,2	0,81	0,96	82,5	1,30
114	0,62	0,55	111	1,68

**Sensitivity:** 1 mmol/L = 0,00471 A.

**Accuracy:** Results obtained using SPINREACT reagents did not show systematic differences when compared with other commercial reagents. The results obtained using 50 samples were the following:

Correlation coefficient (r)<sup>2</sup>: 0,96731.

Regression equation: y=0,990x + 0,100.

The results of the performance characteristics depend on the analyzer used.

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