Quantitative determination of Ceruloplasmin IVD

Store at 2-8°C.

INTENDED USE

The Ceruloplasmin reagent is a quantitative turbidimetric test for the measurement of Ceruloplasmin in human serum or plasma.

PRINCIPLE OF THE METHOD

Anti-human Ceruloplasmin antibodies when mixed with samples containing Ceruloplasmin, form insoluble complexes. These complexes cause an absorbance change, dependent upon the Ceruloplasmin concentration of the patient sample, that can be quantified by comparison from a calibrator of known Ceruloplasmin concentration.

CLINICAL SIGNIFICANCE

Caeruloplasmin is an α_2 -globulin that contains approximately 95% of total serum copper. Each molecule of caeruloplasmin contain six to eight copper atoms. The high content of copper ions gives caeruloplasmin a blue color. Caeruloplasmin can also bind, and possible transport, other cations such as magnesium. The molecule of caeruloplasmin has a single polypeptide chain and carbohydrate, and results a molecular mass of 132 kD. Caeruloplasmin is synthesized primarily by the hepatic cells and small quantities by macrophages and lymphocytes.

Caeruloplasmin is most often quantified as a screening test for Wilson's disease. However, it is important to realize that several other factors, including diet, hormone levels, and other genetic disorders, can influence plasma levels. Synthesis of caeruloplasmin is increased modestly in the acute-phase response, peaking at 4 to 20 days after a single, cute insult. Synthesis is also stimulated by estrogens, and during pregnancy.

Low plasma caeruloplasmin levels are due to lack of incorporation of Cu^{2+} into the molecule during synthesis. The causes are, the dietary insufficiency (including malabsorption), inability to release Cu^{2+} from gastrointestinal epithelium into circulation, or inability to insert Cu^{2+} into developing caeruloplasmin molecule. Levels may also be low in blood loss or gastrointestinal or renal proteinlosing syndromes.

REAGENTS

Diluent (R1)	Tris buffer 20 mmol/L, PEG 8000, pH 8,3. Preservative.
Antibody (R2)	Goat serum, anti-human Ceruloplasmin, pH 7,5. Preservative.
Optional	Cod: 1102003 PROT CAL.

CALIBRATION

The assay is calibrated to the Reference Material CRM 470/RPPHS (Institute for Reference Materials and Measurements). It must be used the PROT CAL to calibrate the reagent.

STORAGE AND STABILITY

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C and contaminations are prevented during their use. Do not use reagents over the expiration date.

Reagent deterioration: The presence of particles and turbidity. Do not use.

Do not freeze; frozen Antibody or Diluent could change the functionality of the test.

ADDITIONAL EQUIPMENT

- Thermostatic bath at 37°C.
- SPIN 800 Autoanalyzer.

SAMPLES

Fresh serum or plasma. EDTA or heparin should be used as anticoagulant. Stable 7 days at 2-8°C or 3 months at -20°C.

The samples with presence of fibrin should be centrifuged.

Do not use highly hemolized or lipemic samples.

REFERENCE VALUES

Between 15 – 60 mg/dL. Each laboratory should establish its own reference range.

QUALITY CONTROL

Control sera are recommended to monitor the performance of manual and automated assay procedures. Spinreact PROT CONTROL (Cod:1102004). Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

BARCODED REAGENTS LOAD MUST BE PRECEDED OF A SPINREACT "DATABASE" COPY INTO THE ANALYZER SOFTWARE. IT IS AVAILABLE UNDER REQUEST TO SPINREACT.

SPIN 800 APPLICATION

EDIT PARAMETERS

Test	CERULOPLASMIN	No.	**
Full name	CER	Print name	CER
Reac. Type	End Point	Direction	Increase
Pri. Wave.	340	Sec. Wave.	
Unit	mg/dL	Decimal	0.1
Reagent Blank	15 -16	React. Time	24- 25
Vol. Sample	2.1 ul	Increases.	Decrease
			Sample blank
R1	240 ul	R2	60 ul
		R3	R4

CALIBRATION

Spline	Water	0	Vol aspirated	Vol net	Vol dil
Calibrator	*0.1		2.1	10	90
Calibrator	*0.25		2.1	30	90
Calibrator	*0.496		4.2	30	90
Calibrator	*0.735		6.3	30	90
Calibrator	(*)Cal value		2.1	--	--

CAL AUTO DILUTIONS

The Calibration is stable until **7 days**. After this period the Calibration must be performed again in order to obtain good results.

PERFORMANCE CHARACTERISTICS

- Linearity limit:** Up to 91 mg/dL, under the described assay conditions. Samples with higher concentrations, should be diluted 1/5 in NaCl 9 g/L and retested again. The linearity limit depends on the sample / reagent ratio. It will be higher by decreasing the sample volume, although the sensitivity of the test will be proportionally decreased.
- Detection Limit:** Values less than 1,12 mg/dL give non-reproducible results.
- Precision:** The reagent has been tested for 20 days, using three levels of serum in a EP5-based study (NCCLS).

EP5	CV (%)		
	28,96mg/dL	55,47 mg/dL	76,54 mg/dL
Total	4%	2,3%	2%
Within Run	2,2%	1,5%	1%
Between Run	3,1%	1,1%	1,5%
Between Day	1,1%	1,3%	0,8%

- Accuracy:** Results obtained using this reagent (y) were compared to those obtained with a Bayer immunoturbidimetric method. 45 samples ranging from 20 to 80 mg/dL of Ceruloplasmin were assayed. The correlation coefficient (r) was 0,96 and the regression equation $y = 0,896x + 10,57$.

The results of the performance characteristics depend on the used analyzer.

BIBLIOGRAPHY

- Clinical Guide to Laboratory Tests, Edited by NW Tietz W B Saunders Co., Philadelphia, 483, 1983.
- Pesce AJ and Kaplan, LA. Methods in Clinical Chemistry. The CV Mosby Company, St. Louis MO, 1987.
- Dati F et al. Eur J Clin Chem Clin Biochem 1966; 14: 401-406.
- Young DS. Effects of disease on clinical laboratory tests, 3th ed. AACC Pres, 1997
- Friedman and Young. Effects of disease on clinical laboratory tests, 3th ed. AACC Pres, 1997.

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BIBLIOGRAPHIE

1. Clinical Guide to Laboratory Tests, Edited by NW Tietz W B Saunders Co., Philadelphia, 483, 1983.
2. Pesce AJ and Kaplan, LA. Methods in Clinical Chemistry. The CV Mosby Company, St. Louis MO, 1987.
3. Dati F et al. Eur J Clin Chem Clin Biochem 1966; 14: 401-406.
4. Young DS. Effects of disease on clinical laboratory tests, 3th ed. AACC Pres, 1997.
5. Friedman and Young. Effects of disease on clinical laboratory tests, 3th ed. AACC Pres, 1997.

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