### Quantitative determination

#### Direct Method

Store at room temperature.  
For "in vitro" diagnostic use only.

#### PRINCIPLE

This method has been developed for the counting of platelets using normal laboratory equipment and whole blood samples.

The technique involves a quick lysis of red cells by the improved diluting in which an effective antiaggregating agent has been added to avoid the shortcomings of conventional methods.

#### REAGENT PREPARATION AND STABILITY

The reagent is ready to use.

Stable when protected from light and stored at room temperature until the expiration date on the label.

Do not use reagents over the expiration date.

#### SAMPLES

Venous blood collected into K<sub>3</sub> – EDTA

#### MATERIALS

- Thoma counting chamber or similar.
- Microscopy

#### PROCEDURE

Ref: 1700091 50 x 1.98 mL

1. Expell into the plastic disposable tube 20 µL (0.02 mL) of blood cells. Mix well by repeated inversion manually or mechanically for 4-5 minutes. Do not shake.
2. Fill the chamber of Thoma by placing the tip of pipette are one of the open ends of the chamber. If air bubbles are present, clean again and refill the chamber.
3. Allow to stand for 15 minutes in order to ensure sedimentation of the platelets in closed Petri dish kept humid.
4. **Counting.** Count total number of platelets of the whole chamber (1 mm<sup>2</sup>). Platelets will appear as small refractile bodies under the X40 objective and the X10 eyepiece.

Ref: 1700090 50 mL

1. Pour into a disposable test tube enough volume to perform the test needed (1mL / Test). Never introduce the red cell pipette directly into a Platelet Counting Fluid bottle.
2. Fill a red cell pipette rapidly with whole blood taken in the way described, to the 1 mark.

3. Diluting fluid is quickly drawn to the 101 mark (dilution 1:100).
4. Shake the pipette either manually or mechanically for 4-5 minutes. Do not shake.
5. Discard the first 3 or 4 drops and then fill the chamber of Thoma by placing the tip of pipette are one of the open ends of the chamber. If air bubbles are present, clean again and refill the chamber.
6. Allow to stand for 15 minutes in order to ensure sedimentation of the platelets in closed Petri dish kept humid.
7. **Counting.** Count total number of platelets of the whole chamber (1 mm<sup>2</sup>). Platelets will appear as small refractile bodies under the X40 objective and the X10 eyepiece.

#### CALCULATIONS

Number of cells counted x 1000 = platelets/mm<sup>3</sup> blood.

#### NORMAL VALUES

150.000 to 400.000 platelet / mm<sup>3</sup>

150 to 400 x 10<sup>9</sup> / L

#### NOTES

The direct counting method should always be supplemented by a blood smear, first, to confirm the counted number and, second, to study platelet morphology.

All thrombocytopenic results must be confirmed before being reported.

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