SPINTROL H NORMAL

SPINTROL H PATHOLOGIC

Human Control serum

SPINTROL H NORMAL

SPINTROL H PATOLÓGICO

Suero Control humano

**Serum control multicomponent for clinical chemistry assays
IVD**

Store at 2 - 8°C.

PRODUCT CHARACTERISTICS

SPINTROL H NORMAL or PATHOLOGIC is a human lyophilised serum. With most constituent concentrations and activities in the normal or pathologic range. It's intended for control of accuracy for use with manual and automated analytical procedures.

REAGENTS

Human serum. Biological additives. Bacteriostatics agents. The concentration/ activities of the components are lot-specific. The exact values and ranges valid for reagents are given informational purpose in the value sheet.

PRECAUTIONS

The individual units comprising lot or lots of components from human origin, have been tested for HBsAg and antibodies to HIV and HCV found negative / non reactive by an FDA approved test. However, handle cautiously as potentially infectious.

PREPARATION

-Reconstitute (→) with 5,0 mL of distilled water.
-Mix thoroughly, avoiding foam forming.
-Bring to room temperature for about 30 min. before use.
Improper handling and/or storage can affect results. Inaccurate reconstitution and errors in assay technique can cause erroneous results.

STORAGE AND STABILITY

The Control serum is stable until the expiration date on the label when stored tightly closed at 2-8°C and contaminations are prevented during their use. Do not use reagents over the expiration date or if there is visible evidence of microbial growth.

Store tightly capped when not use.

-After reconstitution is stable for:

At 15°C to 25°C 12 hours. At 2°C to 8°C 5 days. At -25°C to -15°C 28 days (when frozen once).

-Bilirubin and Acid Phosphatase (stored protected from light):

At 15°C to 25°C 8 hours. At 2°C to 8°C 24 hours. At -25°C to -15°C 2 weeks (when frozen once).

-ALT:

At 15°C to 25°C 12 hours. At 2°C to 8°C 5 days. At -25°C to -15°C 2 weeks (when frozen once).

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