TBA / CO₂ CONTROL

Total Bile Acids / Carbon Dioxide Control-
Control de Ácidos Biliares Totales / Dióxido de
Carbono

Total Bile Acids/ Carbon Dioxide (TBA / CO₂) Control serum

IVD
Store at 2 - 8°C

PRODUCT CHARACTERISTICS

The TBA / CO₂ Control is a bovine lyophilised serum. It is used for quality control procedures in examining the accuracy of quantitative TBA / CO₂ assays.

REAGENTS

Bovine Serum. Biological additives. Bacteriostatics agents. The concentration / activities of the components are lot-specific. The values for each parameter are given on the below table.

PREPARATION

-Reconstitute (→) with 5.0 mL of distilled water.
-Mix thoroughly, avoiding foam forming.
-Bring to room temperature for about 30 min. before use.
Improper handling and/or storage can affect results. Inaccurate reconstitution and errors in assay technique can cause erroneous results.

STORAGE AND STABILITY

The control is stable until the expiration date on the label when stored tightly closed at 2-8°C and contaminations are prevented during their use. Do not use reagents over the expiration date or if there is visible evidence of microbial growth. Store tightly capped when not use.

After reconstitution is stable for:

- At 15°C to 25°C 24 hours
- At 2°C to 8°C 7 days
- At -25°C to -15°C 1 month

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*This Values sheet is applicable to lot and sublots. Sequential alphabetical letter (e.g. A, B, C etc.) following the lot number.

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