

**Quantitative determination of aspartate aminotransferase GOT (AST) IVD**

Store at 2-8°C

**PRINCIPLE OF THE METHOD**

Aspartate aminotransferase (AST) formerly called glutamate oxaloacetate (GOT) catalyses the reversible transfer of an amino group from aspartate to  $\alpha$ -ketoglutarate forming glutamate and oxalacetate. The oxalacetate produced is reduced to malate by malate dehydrogenase (MDH) and NADH:

The rate of decrease in concentration of NADH, measured photometrically, is proportional to the catalytic concentration of AST present in the sample<sup>1</sup>.

**CLINICAL SIGNIFICANCE**

The AST is a cellular enzyme, is found in highest concentration in heart muscle, the cells of the liver, the cells of the skeletal muscle and in smaller amounts in other weaves.

Although an elevated level of AST in the serum is not specific of the hepatic disease, is used mainly to diagnostic and to verify the course of this disease with other enzymes like ALT and ALP.

Also, it is used to control the patients after myocardial infarction, in skeletal muscle disease and other<sup>1,4,5</sup>.

Clinical diagnosis should not be made on a single test result; it should integrate clinical and other laboratory data.

**REAGENTS**

<b>R 1</b>	TRIS pH 7.8	80 mmol/L
<b>Buffer</b>	Lactate dehydrogenase (LDH)	800 U/L
	Malate dehydrogenase (MDH)	600 U/L
	L-Aspartate	200 mmol/L
	<b>R 2</b>	NADH
<b>Substrate</b>	$\alpha$ -Ketoglutarate	12 mmol/L

**PRECAUTIONS**

R1: H290-May be corrosive to metals.

Follow the precautionary statements given in MSDS and label of the product.

**PREPARATION**

All the reagents are ready to use.

**STORAGE AND STABILITY**

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C, protected from light and contaminations prevented during their use.

Do not use reagents over the expiration date.

**Signs of reagent deterioration:**

- Presence of particles and turbidity.
- Blank absorbance (A) at 340 nm < 1,00.

**ADDITIONAL EQUIPMENT**

- Autoanalyzer Spintech 240.
- General laboratory equipment.

**SAMPLES**

Serum or plasma<sup>1</sup>: Stability 7 days at 2-8°C.

**REFERENCE VALUES<sup>1</sup>**

	25°C	30°C	37°C
Men	up to 19 U/L	26 U/L	38 U/L
Women	up to 16 U/L	22 U/L	31 U/L

These values are for orientation purpose; each laboratory should establish its own reference range.

**Temperature conversion factors**

To correct results to other temperatures multiply by:

Assay temperature	Conversion factor to		
	25°C	30°C	37°C
25°C	1,00	1,37	2,08
30°C	0,73	1,00	1,54
37°C	0,48	0,65	1,00

**APPLICATION SPINTECH 240**

Item Name GOT		CALIBRATION	
DATA INFORMATION		TYPE	Linear
Units	U/L		
Decimals	0		
ANALYSIS		STANDARD	
Type	RATE	#1 *	#4
W.Length 1	340	#2	#5
W.Length 2		#3	#6
Method KINETIC		NORMAL RANGE (37°C)	
CORR SLOPE		SERUM MALE	LOW HIGH
1.000 x + 0		URINE FEMALE	
Item Name GOT		DATA PROCESS	
ASPIRATION		ABSORBANCE LIMIT	
KIND	Single <input checked="" type="checkbox"/> Double	READ	LOW -3,000
		START END	HIGH 3,000
SAMPLE	VOLUME	MAIN 37 49	
REAGENT 1	25 $\mu$ L	SUB	
REAGENT 2	200 $\mu$ L		ENDPOINT LIMIT
	50 $\mu$ L		LINEAR CHECK (%) 90
Third Mix	<input checked="" type="checkbox"/> OFF <input type="checkbox"/> ON	FACTOR	
R1 Blank	<input checked="" type="checkbox"/> Water <input type="checkbox"/> R1-B	Blank Correction	1,000
MONITOR		PROZONE CHECK	
0 LEVEL POINT	1	START END	LIMIT (%)
SPAN	3,000	FIRST	
		SECOND	<input checked="" type="checkbox"/> Low <input type="checkbox"/> High
		THIRD	<input checked="" type="checkbox"/> Low <input type="checkbox"/> High

Blank parameter must be performed in order to get good results in CALIB screen from main menu.

**QUALITY CONTROL**

Control sera are recommended to monitor the performance of assay procedures: SPINCONTROL H Normal and Pathologic (Ref. 1002120 and 1002210).

If control values are found outside the defined range, check the instrument, reagents and technique for problems.

Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

**PERFORMANCE CHARACTERISTICS**

**Measuring range:** From detection limit of 0 U/L to linearity limit of 467 U/L.

If the results obtained were greater than linearity limit, dilute the sample 1/10 with NaCl 9 g/L and multiply the result by 10.

**Precision:**

	Intra-assay (n=20)		Inter-assay (n=20)	
	Mean (U/L)	SD	CV (%)	
Mean (U/L)	48,1	159	47,4	156
SD	0,56	0,57	1,42	4,35
CV (%)	1,16	0,36	3,00	2,79

**Sensitivity:** 1 U/L = 0,00053  $\Delta$ A/min.

**Accuracy:** Results obtained using SPINREACT reagents (y) did not show systematic differences when compared with other commercial reagents (x).

The results obtained using 50 samples were the following:

Correlation coefficient (r): 0,99956.

Regression equation: y= 1,042x - 0,342.

The results of the performance characteristics depend on the analyzer used.

**BIBLIOGRAPHY**

- Murray R. Aspartate aminotransferase. Kaplan A et al. Clin Chem The C.V. Mosby Co. St Louis. Toronto. Princeton 1984; 1112-1116.
- Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.
- Young DS. Effects of disease on Clinical Lab. Tests, 4th ed AACC 2001.
- Burtis A et al. Tietz Textbook of Clinical Chemistry, 3rd ed AACC 1999.
- Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed AACC 1995

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**APPLICATION AU SPINTECH 240**

Item Name GOT			
<u>DATA INFORMATION</u>		<u>CALIBRATION</u>	
Units	U/L	TYPE	Linear
Decimals	0		
<u>ANALYSIS</u>		STANDARD	
Type	RATE	#1	*
		#2	#4
W.Length 1	340	#3	#5
W.Length 2			#6
Method KINETIC		<u>NORMAL RANGE (37°C)</u>	
		SERUM	LOW HIGH
<u>CORR</u>		MALE	
SLOPE	INTER	FEMALE	
1.000 x +	0	URINE	
Item Name GOT			
<u>ASPIRATION</u>		<u>DATA PROCESS</u>	
KIND	Single <input type="checkbox"/> Double <input checked="" type="checkbox"/>	<u>READ</u>	LOW -3,000
		START END	HIGH 3,000
SAMPLE	VOLUME	MAIN 37 49	
REAGENT 1	25 µL	SUB	
REAGENT 2	200 µL		
REAGENT 2	50 µL		ENDPOINT LIMIT
			LINEAR CHECK (%) 90
Third Mix	<input checked="" type="checkbox"/> OFF <input type="checkbox"/> ON	<u>FACTOR</u>	
R1 Blank	<input checked="" type="checkbox"/> Water <input type="checkbox"/> R1-B	Blank Correction	1,000
<u>MONITOR</u>		<u>PROZONE CHECK</u>	
0 LEVEL POINT	1	START END	LIMIT (%)
SPAN	3,000	FIRST	<input type="checkbox"/> Low <input type="checkbox"/> High
		SECOND	<input type="checkbox"/> Low <input type="checkbox"/> High
		THIRD	<input type="checkbox"/> Low <input type="checkbox"/> High

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Item Name GOT <b>DATA INFORMATION</b> Units U/L Decimals 0 <b>ANALYSIS</b> Type RATE W.Length 1 340 W.Length 2 Method KINETIC <b>CORR</b> SLOPE INTER 1.000 x + 0		<b>CALIBRATION</b> TYPE Linear STANDARD #1 * #4 #2 #5 #3 #6 <b>NORMAL RANGE (37°C)</b> LOW HIGH SERUM MALE FEMALE URINE	
Item Name GOT <b>ASPIRATION</b> KIND Single <input checked="" type="checkbox"/> Double VOLUME SAMPLE 25 µL REAGENT 1 200 µL REAGENT 2 50 µL Third Mix <input checked="" type="checkbox"/> OFF ON R1 Blank <input checked="" type="checkbox"/> Water R1-B		<b>DATA PROCESS</b> <b>READ</b> START END MAIN 37 49 SUB ENDPOINT LIMIT LINEAR CHECK (%) 90 <b>FACTOR</b> Blank Correction 1,000 <b>PROZONE CHECK</b> START END LIMIT (%) FIRST SECOND <input checked="" type="checkbox"/> Low High THIRD <input checked="" type="checkbox"/> Low High	
<b>MONITOR</b> 0 LEVEL POINT 1 SPAN 3,000			

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