
**Quantitative determination of gamma-glutamyl transferase (γ-GT) IVD**

Store at 2-8°C

**PRINCIPLE OF THE METHOD**

Gamma-glutamyl transferase (γ-GT) catalyses the transfer of γ-glutamyl group from γ-glutamyl-p-nitroanilide to acceptor glycylglycine, according to the following reaction:

The rate of 2-nitro-5-aminobenzoic acid formation, measured photometrically, is proportional to the catalytic concentration of γ-GT present in the sample<sup>1,2</sup>.

**CLINICAL SIGNIFICANCE**

Gamma-glutamyl transferase (γ-GT) is a cellular enzyme with wide tissue distribution in the body, primarily in the kidney, pancreas, liver and prostate. Measurements of gamma-glutamyl transferase (γ-GT) activity are used in the diagnosis and treatment of hepatobiliary diseases such as biliary obstruction, cirrhosis or liver tumours<sup>1,2,5,6</sup>.

Clinical diagnosis should not be made on a single test result; it should integrate clinical and other laboratory data.

**REAGENTS**

<b>R 1</b>	TRIS	100 mmol/L
Buffer	Glycylglycine	100 mmol/L
<b>R 2</b>	L-γ-glutamyl-3-carboxy-4-nitroanilide	3 mmol/L
Substrate		

**PREPARATION**

All the reagents are ready to use.

**STORAGE AND STABILITY**

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C, protected from light and contaminations prevented during their use.

Do not use reagents over the expiration date.

**Signs of reagent deterioration:**

- Presence of particles and turbidity.
- Blank absorbance (A) at 405 nm ≥ 1,20.

**ADDITIONAL EQUIPMENT**

- Autoanalyzer Spintech 240.
- General laboratory equipment.

**SAMPLES**

Serum<sup>1</sup>. γ-GT is stable for at least 3 days at 2-8°C, 8 hours at 15-25°C and 1 month at -20°C.

**REFERENCE VALUES<sup>1</sup>**

	25°C	30°C	37°C
Women	4-18 U/L	5-25 U/L	7-32 U/L
Men	6-28 U/L	8-38 U/L	11-50 U/L

These values are for orientation purpose; each laboratory should establish its own reference range.

**Temperature conversion factors**

To correct results to other temperatures multiply by:

Assay temperature	Conversion factor to		
	25°C	30°C	37°C
25°C	1,00	1,37	1,79
30°C	0,73	1,00	1,30
37°C	0,56	0,77	1,00

**APPLICATION SPINTECH 240**

<p>Item Name GGT</p> <p><u>DATA INFORMATION</u></p> <p>Units U/L Decimals 0</p> <p><u>ANALYSIS</u></p> <p>Type RATE</p> <p>W.Length 1 405 W.Length 2 700</p> <p>Method Carboxy</p> <p><u>CORR</u></p> <p>SLOPE INTER 0 1.000 x +</p>	<p><u>CALIBRATION</u></p> <p>TYPE Linear</p> <p>STANDARD #1 * #4 #2 #5 #3 #6</p> <p><u>NORMAL RANGE (37°C)</u></p> <p>SERUM MALE LOW HIGH FEMALE</p> <p>URINE</p>
<p>Item Name GGT</p> <p><u>ASPIRATION</u></p> <p>KIND Single <input checked="" type="checkbox"/> Double</p> <p>VOLUME</p> <p>SAMPLE 25 µL REAGENT 1 200 µL REAGENT 2 50 µL</p> <p>Third Mix <input checked="" type="checkbox"/> OFF ON R1 Blank <input checked="" type="checkbox"/> Water R1-B</p> <p><u>MONITOR</u></p> <p>0 LEVEL POINT 1 SPAN 3.000</p>	<p><u>DATA PROCESS</u></p> <p><u>READ</u></p> <p>START END LOW HIGH MAIN 36 44 -3,000 3,000 SUB</p> <p><u>ENDPOINT LIMIT</u></p> <p>LINEAR CHECK (%) 90</p> <p><u>FACTOR</u></p> <p>Blank Correction 1,000</p> <p><u>PROZONE CHECK</u></p> <p>START END LIMIT (%)</p> <p>FIRST SECOND <input checked="" type="checkbox"/> Low High THIRD <input checked="" type="checkbox"/> Low High</p>

*Blank parameter must be performed in order to get good results in CALIB screen from main menu.*

**QUALITY CONTROL**

Control sera are recommended to monitor the performance of assay procedures: SPINTROL H Normal and Pathologic (Ref. 1002120 and 1002210).

If control values are found outside the defined range, check the instrument, reagents and technique for problems.

Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

**PERFORMANCE CHARACTERISTICS**

**Measuring range:** From detection limit of 2 U/L to linearity limit of 300 U/L.

If the results obtained were greater than linearity limit, dilute the sample 1/10 with NaCl 9 g/L and multiply the result by 10.

**Precision:**

	Intra-assay (n=20)		Inter-assay (n=20)	
	Mean (U/L)	SD	Mean (U/L)	SD
Mean (U/L)	38,3	190	40,1	198
SD	0,39	0,53	0,82	2,30
CV (%)	1,03	0,28	2,05	1,16

**Sensitivity:** 1 U/L = 0,0008 ΔA/min.

**Accuracy:** Results obtained using SPINREACT reagents (y) did not show systematic differences when compared with other commercial reagents (x).

The results obtained using 50 samples were the following:

Correlation coefficient (r)<sup>2</sup>: 0,99990.

Regression equation: y= 1,334x - 1,493.

The results of the performance characteristics depend on the analyzer used.

**BIBLIOGRAPHY**

1. Gendler S. γ-GT. Kaplan A et al. Clin Chem The C.V. Mosby Co. St Louis. Toronto. Princeton 1984; 1120-1123.
2. Persijn J P et al. J Clin Chem Clin Biochem 1976; (14) 9: 421-427.
3. Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.
4. Young DS. Effects of disease on Clinical Lab. Tests, 4th ed AACC 2001.
5. Burtis A et al. Tietz Textbook of Clinical Chemistry, 3rd ed AACC 1999.
6. Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed AACC 1995.

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