Quantitative determination of albumin
IVD

Store at 2-8°C

PRINCIPLE OF THE METHOD

Albumin in the presence of bromocresol green at a slightly acid pH, produces a colour change of the indicator from yellow-green to green-blue. The intensity of the colour formed is proportional to the albumin concentration in the sample^{1,2,3,4}.

CLINICAL SIGNIFICANCE

One of the most important serum proteins produced in the liver is albumin. This molecule has an extraordinarily wide range of functions, including nutrition, maintenance of oncotic pressure and transport of Ca⁺⁺, bilirubin, free fatty acid, drugs and steroids.

Variation in albumin levels indicate liver diseases, malnutrition, skin lesions such as dermatitis and burns or dehydration^{1,7,8}.

Clinical diagnosis should not be made on a single test result; it should integrate clinical and other laboratory data.

REAGENTS

R	Bromocresol green pH 4.2	0.12 mmol/L
----------	--------------------------	-------------

PREPARATION

The reagent is ready to use.

STORAGE AND STABILITY

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C, protected from light and contaminations prevented during their use.

Do not use reagents over the expiration date.

Signs of reagent deterioration:

- Presence of particles and turbidity.
- Blank absorbance (A) at 600 nm \geq 0.40.

ADDITIONAL EQUIPMENT

- Autoanalyzer Spintech 240.
- General laboratory equipment.

SAMPLES

Serum or plasma, free of hemolysis¹: Stability 1 month at 2-8°C or 1 week at 15-25°C.

REFERENCE VALUES

 3.5 to 5.0 g/dL¹.

These values are for orientation purpose; each laboratory should establish its own reference range.

Conversion factor: g/dL x 144.9 = μ mol/L

QUALITY CONTROL

Control sera and calibrators are recommended to monitor the performance of assay procedures: SPINTROL H Calibrator, SPINTROL H Normal and Pathologic (Ref. 1002011, 1002120 and 1002210).

If control values are found outside the defined range, check the instrument, reagents and technique for problems.

Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

NOTES

1. Calibration with the aqueous Standard may cause a systematic error in automatic procedures. In these cases, it is recommended to use a serum Calibrator.
2. Use clean disposable pipette tips for its dispensation.

APPLICATION SPINTECH 240

Item Name ALB DATA INFORMATION Units g/dL Decimals 2 ANALYSIS Type END W.Length 1 600 Method BCG CORR SLOPE INTER 1.000 x + 0		CALIBRATION TYPE Linear STANDARD #1 * #4 #2 #5 #3 #6 NORMAL RANGE LOW HIGH SERUM MALE FEMALE	
Item Name ALB ASPIRATION KIND <input checked="" type="checkbox"/> Single <input type="checkbox"/> Double VOLUME SAMPLE 3 μ L REAGENT 1 300 μ L Third Mix <input checked="" type="checkbox"/> OFF <input type="checkbox"/> ON R1 Blank Water <input checked="" type="checkbox"/> R1-B		DATA PROCESS ABSORBANCE LIMIT READ START END LOW HIGH MAIN 50 52 -3.000 3.000 SUB ENDPOINT LIMIT 3 LINEAR CHECK (%) FACTOR Blank Correction 1.000 PROZONE CHECK START END LIMIT (%) FIRST SECOND <input checked="" type="checkbox"/> Low High THIRD <input checked="" type="checkbox"/> Low High	
MONITOR 0 LEVEL POINT 1 SPAN 3.000			

Blank parameter must be performed in order to get good results in CALIB screen from main menu. This parameter calibration is stable for more than 40 days

PERFORMANCE CHARACTERISTICS

Measuring range: From detection limit of 0,0349 g/dL to linearity limit of 6 g/dL. If the results obtained were greater than linearity limit, dilute the sample 1/2 with NaCl 9 g/L and multiply the result by 2.

Precision:

	Intra-assay (n=20)		Inter-assay (n=20)	
	Mean (g/dL)	SD	Mean (g/dL)	SD
Mean (g/dL)	5,00	3,71	4,56	3,07
SD	0,02	0,02	0,28	0,18
CV (%)	0,47	0,55	6,20	5,90

Sensitivity: 1 g/dL = 0,2003 A.

Accuracy: Results obtained using SPINREACT reagents (y) did not show systematic differences when compared with other commercial reagents (x).

The results obtained using 50 samples were the following:

Correlation coefficient (r)²: 0,99169.

Regression equation: y= 1,045x – 0,028.

The results of the performance characteristics depend on the analyzer used.

BIBLIOGRAPHY

1. Gendler S. Uric acid. Kaplan A et al. Clin Chem The C.V. Mosby Co. St Louis. Toronto. Princeton 1984; 1268-1273 and 425.
2. Rodkey F L. Clin Chem 1965; 11: 478-487.
3. Webster D. Clin Chem. 1974; Acta 53: 109-115.
4. Doumas BT Clin Chem. 1971; Acta 31: 87-96.
5. Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.
6. Young DS. Effects of disease on Clinical Lab. Tests, 4th ed AACC 2001.
7. Burtis A et al. Tietz Textbook of Clinical Chemistry, 3rd ed AACC 1999.
8. Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed AACC 1995.

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