Quantitative determination of calcium IVD

Store at 2-8°C

PRINCIPLE OF THE METHOD

Calcium with Arsenazo III (1,8-Dihydroxy-3,6-disulpho-2,7-naphthalene-bis (azo)-dibenzeneazonic acid), at neutral pH, yields a blue colored complex.

The intensity of the colour formed is proportional to the calcium concentration in the sample^{1,2,3}.

CLINICAL SIGNIFICANCE

Calcium is the most abundant and one of the most important minerals in the human body. Approximately 99% of body calcium is found in bones.

A decrease in albumin level causes a decrease in serum calcium. Low levels of calcium are found in hypoparathyroidism, pseudohypoparathyroidism, vitamin D deficiency, malnutrition and intestinal malabsorption.

Among causes of hypercalcemia are cancers, large intake of vitamin D, enhanced renal retention, osteoporosis, sarcoidosis, thyrotoxicosis, hyperparathyroidism^{1,6,7}.

Clinical diagnosis should not be made on a single test result; it should integrate clinical and other laboratory data.

REAGENTS

R	Imidazol Buffer pH 6.5	100 mmol/L
Arsenazo III	Arsenazo III	120 mmol/L

PRECAUTIONS

R: H360- May damage fertility or the unborn child.

Follow the precautionary statements given in MSDS and label of the product.

PREPARATION

The reagent is ready to use.

STORAGE AND STABILITY

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C, protected from light and contaminations are prevented during their use.

Do not use reagents over the expiration date.

Signs of reagent deterioration:

- Presence of particles and turbidity.
- Blank absorbance (A) at 660 nm \geq 0.50.

ADDITIONAL EQUIPMENT

- Autoanalyzer Spintech 240.
- General laboratory equipment (Note 1, 2).

SAMPLES

- Serum or plasma¹: Separated from cells as rapidly as possible. Blood anticoagulants with oxalate or EDTA are not acceptable since these chemicals will strongly chelate calcium.

- Urine¹: Collect 24 hours urine specimen in calcium free containers. The collecting bottles should contain 10 ml of diluted Nitric acid (50% v/v). Record the volume.

Dilute a sample 1/2 in distilled water. Mix. Multiply results by 2 (dilution factor).

Stability of the samples: Calcium is stable 10 days at 2-8°C.

REFERENCE VALUES¹

Serum or plasma:

Adults	8.5-10.5 mg/dL	\cong 2.1-2.6 mmol/L
Children	10 -12 mg/dL	\cong 2.5 - 3 mmol/L
Newborns	8 -13 mg/dL	\cong 2 - 3.25 mmol/L

Urine:

Adults	50 - 300 mg/24h	\cong 1.25 - 7.5 mmol/24h
Children	80 -160 mg/24h	\cong 2 - 4 mmol/24h

These values are for orientation purpose; each laboratory should establish its own reference range.

Conversion factor: mg/dL x 0.25= mmol/L.

QUALITY CONTROL

Control sera and calibrators are recommended to monitor the performance of assay procedures: SPINROL H Calibrator, SPINROL H Normal and Pathologic (Ref. 1002011, 1002120 and 1002210).

If control values are found outside the defined range, check the instrument, reagents and technique for problems.

Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

APPLICATION SPINTECH 240

Item Name CA			
<u>DATA INFORMATION</u>			
Units	mg/dL		
Decimals	2		
<u>ANALYSIS</u>			
Type	END		
W.Length 1	660		
Method	Arsenazo III		
<u>CORR</u>			
SLOPE	INTER		
1.000 x +	0		
Item Name CA			
<u>ASPIRATION</u>			
KIND	<input checked="" type="checkbox"/> Single	<input type="checkbox"/> Double	
<u>VOLUME</u>			
SAMPLE	3 μ L		
REAGENT 1	300 μ L		
Third Mix	<input checked="" type="checkbox"/> OFF	<input type="checkbox"/> ON	
R1 Blank	Water	<input checked="" type="checkbox"/> R1-B	
<u>MONITOR</u>			
0 LEVEL POINT	1		
SPAN	3.000		
		<u>CALIBRATION</u>	
		TYPE	Linear
		<u>STANDARD</u>	
		#1	* #4
		#2	#5
		#3	#6
		<u>NORMAL RANGE</u>	
		LOW	HIGH
		SERUM	MALE
			FEMALE
		<u>DATA PROCESS</u>	
		<u>READ</u>	<u>ABSORBANCE LIMIT</u>
		START	END
		MAIN 30	31
		SUB	
		<u>ENDPOINT LIMIT</u> 3	
		<u>LINEAR CHECK (%)</u>	
		<u>FACTOR</u>	
		Blank Correction	1.000
		<u>PROZONE CHECK</u>	
		START	END
		FIRST	LIMIT (%)
		SECOND	<input checked="" type="checkbox"/> Low High
		THIRD	<input checked="" type="checkbox"/> Low High

Blank parameter must be performed in order to get good results in CALIB screen from main menu. This parameter calibration is stable for more than 40 days.

PERFORMANCE CHARACTERISTICS

Measuring range: From detection limit of 0,026 mg/dL to linearity limit of 32 mg/dL.

If the results obtained were greater than linearity limit, dilute the sample 1/2 with NaCl 9 g/L and multiply the result by 2.

Precision:

	Intra-assay (n=20)	
Mean (mg/dL)	8,35	14,28
SD	0,08	0,08
CV (%)	0,95	0,59

	Inter-assay (n=20)	
	8,58	14,57
	0,19	0,34
	2,24	2,31

Sensitivity: 1 mg/dL = 0,0316 A.

Accuracy: Results obtained using SPINREACT reagents (y) did not show systematic differences when compared with other commercial reagents (x).

The results obtained using 50 samples were the following:

Correlation coefficient (r): 0,9506

Regression equation: $y = 0,8944x + 1,3421$.

The results of the performance characteristics depend on the analyzer used.

NOTES

1. It is recommended to use disposable material. If glassware is used the material should be scrupulously cleaned with diluted 1/1 HNO₃ in water and then thoroughly rinsed it with distilled water.
2. Most of the detergents and water softening products used in the laboratories contain chelating agents. A defective rinsing will invalidate the procedure.
3. Calibration with the aqueous standard may cause a systematic error in automatic procedures. In these cases, it is recommended to use a serum Calibrator.
4. Use clean disposable pipette tips for its dispensation.

BIBLIOGRAPHY

1. Farell E C. Calcium. Kaplan A et al. Clin Chem The C.V. Mosby Co. St Louis. Toronto. Princeton 1984; 1051-1255 and 418.
2. Kessler G. et al. Clin Chem Vol 10, No 8 1964; 686-706.
3. Connerty H. V. et al. Am J Clin Path Vol 45, No 3 1996; 200-296.
4. Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.
5. Young DS. Effects of disease on Clinical Lab. Tests, 4th ed. AACC 2001.
6. Burtis A. et al. Tietz Textbook of Clinical Chemistry, 3rd ed. AACC 1999.
7. Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed. AACC 1995.

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BIBLIOGRAFIA

1. Farell E C. Calcium. Kaplan A et al. Clin Chem The C.V. Mosby Co. St Louis. Toronto. Princeton 1984; 1051-1255 and 418.
2. Kessler G. et al. Clin Chem 1964; 10 (8); 686-706.
3. Connerty H. V. et al. Am J Clin Path 1996; 45 (3); 200-296.
4. Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.
5. Young DS. Effects of disease on Clinical Lab. Tests, 4th ed. AACC 2001.
6. Burtis A. et al. Tietz Textbook of Clinical Chemistry, 3rd ed. AACC 1999.
7. Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed. AACC 1995.

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BIBLIOGRAPHIE

- Farell E C. Calcium. Kaplan A et al. Clin Chem The C.V. Mosby Co. St Louis. Toronto. Princeton 1984; 1051-1255 and 418.
- Kessler G. et al. Clin Chem 1964; 10 (8); 686-706.
- Connerty H. V. et al. Am J Clin Path 1996; 45 (3); 200-296.
- Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.
- Young DS. Effects of disease on Clinical Lab. Tests, 4th ed. AACC 2001.
- Burtis A. et al. Tietz Textbook of Clinical Chemistry, 3rd ed. AACC 1999.
- Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed. AACC 1995.

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