
Quantitative determination of chloride ion IVD

Store at 2-8°C

PRINCIPLE OF THE METHOD

The quantitative displacement of thiocyanate by chloride from mercuric thiocyanate and subsequent formation of a red ferric thiocyanate complex is measured colorimetrically:

The intensity of the color formed is proportional to the chloride ion concentration in the sample^{1,2,3,4}.

CLINICAL SIGNIFICANCE

It is important clinically the determination of chloride due regulation of osmotic pressure of extra cellular fluid and to its significant role in acid-base balance. Increases in chloride ion concentration may be found in severe dehydration, excessive intake of chloride, severe renal tubular damage and in patients with cystic fibrosis.

Decrease in chloride ion concentration may be found in metabolic acidosis, loss from prolonged vomiting and chronic pyelonephritis^{2,7,8}.

Clinical diagnosis should not be made on a single test result; it should integrate clinical and other laboratory data.

REAGENTS

R Thiocyanate-Hg	Mercuric thiocyanate	4 mmol/L
	Ferric nitrate	40 mmol/L
	Mercuric nitrate	2 mmol/L
	Nitric acid	45 mmol/L

PRECAUTIONS

R: H314-Causes severe skin burns and eye damage.

Follow the precautionary statements given in MSDS and label of the product.

PREPARATION

The reagent is ready to use.

STORAGE AND STABILITY

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C, protected from light and contaminations prevented during their use. Do not use reagents over the expiration date.

Signs of reagent deterioration:

- Presence of particles and turbidity.
- Blank absorbance (A) at 505 nm $\geq 0,15$.

ADDITIONAL EQUIPMENT

- Autoanalyzer Spintech 240.
- General laboratory equipment (Note 1, 2, 3).

SAMPLES

- Serum or plasma free of haemolysis and separated from cells as rapidly as possible. Anticoagulants such as oxalate or EDTA are not acceptable they will interfere with results.

- Urine¹: Collect 24-hour urine specimen in chloride free containers. Dilute a sample 1/2 in distilled water. Mix. Multiply results by 2 (dilution factor).

Stability of the sample: Ion chloride is stable 1 week at room temperature (15-25°C), in refrigerator (2-8°C) or frozen (-20°C) temperatures.

REFERENCE VALUES¹

Serum or plasma:	95 - 115 mmol/L
Urine:	110 - 250 mmol/24h

These values are for orientation purpose; each laboratory should establish its own reference range.

Conversion factor: mmol/L= mEq/L.

QUALITY CONTROL

Control sera and calibrators are recommended to monitor the performance of assay procedures: SPINROL H Calibrator, SPINROL H Normal and Pathologic (Ref. 1002011, 1002120 and 1002210).

If control values are found outside the defined range, check the instrument, reagents and technique for problems.

Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

APPLICATION SPINTECH 240

Item Name Chloride		CALIBRATION	
DATA INFORMATION		TYPE Linear	
Units	mmol/L		
Decimals	1		
ANALYSIS		STANDARD	
Type	END	#1	* #4
		#2	#5
		#3	#6
W.Length 1	505	NORMAL RANGE	
Method Thiocyanate-Hg		LOW HIGH	
CORR SLOPE		SERUM MALE FEMALE	
1.000 x +	INTER 0		
Item Name Chloride		DATA PROCESS	
ASPIRATION		ABSORBANCE LIMIT	
KIND	Single Double	READ	LOW -3.000
		START END	HIGH 3.000
SAMPLE	VOLUME	MAIN 30 31	
REAGENT 1	300 µL	SUB	ENDPOINT LIMIT 3
			LINEAR CHECK (%)
Third Mix		FACTOR	
R1 Blank	OFF ON	Blank Correction	1.000
	Water R1-B	PROZONE CHECK	
MONITOR		START END LIMIT (%)	
0 LEVEL POINT	1	FIRST	
SPAN	3.000	SECOND	Low High
		THIRD	Low High

Blank parameter must be performed in order to get good results in CALIB screen from main menu. This parameter calibration is stable for more than 40 days.

PERFORMANCE CHARACTERISTICS

Measuring range: From detection limit of 0,454 mmol/L to linearity limit of 190 mmol/L.

If the results obtained were greater than linearity limit, dilute the sample 1/2 with distilled water and multiply the result by 2.

Precision:

	Intra-assay (n=20)		Inter-assay (n=20)	
	Mean (mmol/L)	84,2	114	82,5
SD	0,81	0,62	1,07	1,87
CV (%)	0,96	0,55	1,30	1,68

Sensitivity: 1 mmol/L = 0,00471 A.

Accuracy: Results obtained using SPINREACT reagents did not show systematic differences when compared with other commercial reagents. The results obtained using 50 samples were the following:

Correlation coefficient (r)²: 0,96731.

Regression equation: y=0,990x + 0,100.

The results of the performance characteristics depend on the analyzer used.

NOTES

1. It is recommended to use disposable material. If glassware is used the material should be scrupulously cleaned with H₂SO₄ - K₂Cr₂O₇ Solution and then thoroughly rinsed it with distilled water.
2. Avoid the contact with metal materials.
3. Use clean disposable pipette tips for its dispensation.

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