Quantitative determination of Pre-albumin IVD

Store 2 - 8°C.

PRINCIPLE OF THE METHOD

PRE-ALBUMIN is a quantitative turbidimetric test for the measurement of prealbumin in human serum or plasma. Anti-prealbumin antibodies when mixed with samples containing prealbumin, form insoluble complexes. These complexes cause an absorbance change, dependent upon the prealbumin concentration of the patient sample, that can be quantified by comparison from a calibrator of known prealbumin concentration.

CLINICAL SIGNIFICANCE

The prealbumin is a non-glycosylated protein synthesized mainly in the liver and choroid plexus of the brain. It binds and transport approximately 10% of serum thyroxin and triiodothyronine, and also plays a role in the transport of vitamin A in complex with retinal-binding protein.

Prealbumin is the earliest laboratory indicator of nutritional status and has emerged as the preferred marker for malnutrition because it correlates with patient outcomes in wide variety of clinical conditions. It is also a negative acute phase protein; serum levels fall in inflammation and malignancy, as well as cirrhosis, protein-losing enteropathy and zinc deficiency. However, the presence of a prealbumin producing tumor or Hodgkin's disease will increase serum concentrations.

REAGENTS

Diluent (R1)	Tris buffer 20 mmol/L, PEG 8000, pH 8.3. Preservative.
Antibody (R2)	Goat serum, anti-human prealbumin pH 7.5. Preservative.
Optional	Ref: 1102003 PROT CAL

CALIBRATION

The assay has been standardized against the Reference Material ERM-DA470k/IFCC. It is recommended the use of the PROT CAL for calibration.

PREPARATION

Reagents: Ready to use.

Calibration Curve: Prepare the following PROT CAL Calibrator dilutions in NaCl 9 g/L as diluent. Multiply the concentration of the pre-albumin calibrator by the corresponding factor stated in table below to obtain the prealbumin concentration of each dilution.

Calibrator dilution	1	2	3	4	5	6
Calibrator (µL)	--	10	25	50	75	100
NaCl 9 g/L (µL)	100	90	75	50	25	-
Factor	0	0.1	0.25	0.5	0.75	1.0

STORAGE AND STABILITY

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C and contaminations are prevented during their use. Do not use reagents over the expiration date.

Reagent deterioration: The presence of particles and turbidity.

Do not freeze; frozen Antibody or Diluent could change the functionality of the test.

ADDITIONAL EQUIPMENT

- Spintech 240 autoanalyzer
- Laboratory equipment.

SAMPLES

Fresh serum or plasma. EDTA or heparin should be used as anticoagulant. Stable 7 days at 2-8°C or 3 months at -20°C.

The samples with presence of fibrin should be centrifuged before testing.

Do not use highly hemolyzed or lipemic samples.

QUALITY CONTROL

Control sera are recommended to monitor the performance of manual and automated assay procedures. It should be used the SPINREACT PROT CONTROL (Ref.:1102004). Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

REFERENCE VALUES²

Between 20 - 40 mg/dL. Each laboratory should establish its own reference range.

INTERFERENCES

Hemoglobin (16 g/L), bilirubin (40 mg/dL), rheumatoid factors (200 IU/mL), do not interfere. Lipemia (≥ 8 g/L), interfere. Other substances may interfere ^{5,6}.

SPINTECH 240 APPLICATION

Item Name Pre-albumin			
DATA INFORMATION			
Units	mg/dL		
Decimals	0		
ANALYSIS			
Type	END		
W.Length 1	340		
Method	Turbidimetry		
SLOPE	INTER		
1.000 x +	0		
Item Name Pre-albumin			
ASPIRATION			
KIND	Single	<input checked="" type="checkbox"/> Double	
		VOLUME**	
SAMPLE	2	µL	
REAGENT 1	240	µL	
REAGENT 2	60	µL	
Third Mix	<input checked="" type="checkbox"/> OFF	ON	
R1 Blank	<input checked="" type="checkbox"/> Water	R1-B	
MONITOR			
0 LEVEL POINT	1		
SPAN	3.000		
		START	END
FIRST			
SECOND			
THIRD			
		LIMIT (%)	
		<input checked="" type="checkbox"/> Low	<input checked="" type="checkbox"/> High
		<input checked="" type="checkbox"/> Low	<input checked="" type="checkbox"/> High

** Modify reagents and sample volumes according to the range accepted but keeping always the mentioned ratio.

PERFORMANCE CHARACTERISTICS

1. Measurement range: Up to 100 mg/dL, under the described assay conditions. Samples with higher concentrations, should be diluted 1/5 in NaCl 9 g/L and retested again. The linearity limit depends on the sample /reagent ratio. It will be higher by decreasing the sample volume, although the sensitivity of the test will be proportionally decreased.

2. Detection Limit: Values less than 0,69 mg/dL give non-reproducible results.

3. Prozone effect: No prozone effect was detected up to 230 mg/dL.

4. Precision:

EP5	CV (%)		
	20,04mg/dL	39,49 mg/dL	58,67 mg/dL
Total	4,9%	4,6%	4,9%
Within Run	1,8%	2,3%	1,4%
Between Run	2,4%	2,7%	2,4%
Between Day	3,9%	3,0%	4,0%

5. Accuracy: Results obtained using this reagent (y) were compared to those obtained using another immunoturbidimetric method, 60 samples ranging from 1 to 40 mg/dL of pre-albumin were assayed. The correlation coefficient (r) was 0,93 and the regression equation $y = 1,031x - 4,617$,

The results of the performance characteristics depend on the used analyzer,

NOTES

1. Clinical diagnosis should not be made on findings of a single test result, but should integrate both clinical and laboratory data.

BIBLIOGRAPHY

1. Clinical Guide to Laboratory Tests, Edited by NW Tietz W B Saunders Co., Philadelphia, 483, 1983.
2. Dati F et al. Eur J Clin Chem Clin Biochem 1996; 34:517-520.
3. Pesce AJ and Kaplan, LA. Methods in Clinical Chemistry. The CV Mosby Company, St. Louis MO, 1987.
4. Jayle MF et al. Progress in Hematology 1962; 3: 343-359.
5. Young DS. Effects of drugs on clinical laboratory tests, 4th ed. AACC Pres, 1995.
6. Friedman and Young. Effects of disease on clinical laboratory tests, 3th ed. AACC Pres, 1997.

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1. Clinical Guide to Laboratory Tests, Edited by NW Tietz W B Saunders Co., Philadelphia, 483, 1983.
2. Dati F et al. Eur J Clin Chem Biochem 1996; 34:517-520.
3. Pesce AJ and Kaplan, LA. Methods in Clinical Chemistry. The CV Mosby Company, St. Louis MO, 1987.
4. Jayle MF et al. Progress in Hematology 1962; 3: 343-359.
5. Young DS. Effects of drugs on clinical laboratory tests, 4th ed. AACC Pres, 1995.
6. Friedman and Young. Effects of disease on clinical laboratory tests, 3th ed. AACC Pres, 1997.

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BIBLIOGRAPHIE

- Clinical Guide to Laboratory Tests, Edited by NW Tietz W B Saunders Co., Philadelphia, 483, 1983.
- Dati F et al. Eur J Clin Chem Biochem 1996; 34:517-520.
- Pesce AJ and Kaplan, LA. Methods in Clinical Chemistry. The CV Mosby Company, St. Louis MO, 1987.
- Jayle MF et al. Progress in Hematology 1962; 3: 343-359.
- Young DS. Effects of drugs on clinical laboratory tests, 4th ed. AACC Pres, 1995.
- Friedman and Young. Effects of disease on clinical laboratory tests, 3th ed. AACC Pres, 1997.

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BIBLIOGRAFIA

1. Clinical Guide to Laboratory Tests, Edited by NW Tietz W B Saunders Co., Philadelphia, 483, 1983.
2. Dati F et al. Eur J Clin Chem Clin Biochem 1996; 34:517-520.
3. Pesce AJ and Kaplan, LA. Methods in Clinical Chemistry. The CV Mosby Company, St. Louis MO, 1987.
4. Jayle MF et al. Progress in Hematology 1962; 3: 343-359
5. Young DS. Effects of drugs on clinical laboratory tests, 4th ed. AAC Pres, 1995.
6. Friedman and Young. Effects of disease on clinical laboratory tests, 3th ed. AAC Pres, 1997.

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Ref.: TK1102124

Cont.

R1 Diluyente: 2 x 24 mL
R2 Anticorpo: 2 x 6 mL