Quantitative determination of copper IVD

Store at 2-8°C

PRINCIPLE OF THE METHOD

Direct colorimetric test without deproteinization of the sample. End point increase. At pH 4.70, in a buffered media, copper is released from ceruloplasmine complex and forms with the specific complexant 3-5 Di Br-PAESA a stable coloured complex.

The color intensity is proportional to the amount of copper present in the sample^{4,2}.

CLINICAL SIGNIFICANCE

A variety of human copper deficiency conditions are recognized. Specific diseases associated with copper include head disease, bone and joint osteoarthritis and osteoporosis and Menkes' syndrome, Wilson's disease and others. Elevated levels of copper can also be toxic.

Clinical diagnosis should not be made on a single test result; it should integrate clinical and other laboratory data.

REAGENTS

R 1	Acetate. pH 4,7	≥ 1 mol/L
R 2	3,5-DiBr-PAESA	0,4 mmol/L
Color	Acetic acid	1400 mmol/L
	Sodium hydroxide	500 mmol/L
R 3	Ascorbic acid (powder)	
Reducing acid		
COPPER CAL	Copper primary standard	100 µg/dL

PRECAUTIONS

R1: H315-Causes skin irritation. H319-Causes serious eye irritation.

R2: H314-Causes severe skin burns and eye damage.

Follow the precautionary statements given in MSDS and label of the product.

PREPARATION AND STABILITY

- Working reagent (WR): Add (→) one dose (dispense using the enclosed spoon) of R 3 to one vial of R1. Cap and mix gently to dissolve contents. Stability 15 days at 2-8°C when stored tightly closed and contaminations prevented during their use.

- R2: Ready to use. After opening, is stable 90 days 2-8°C, if contamination avoided and vial recapped immediately after use.

- COPPER CAL: Ready to use.

STORAGE AND STABILITY

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C, protected from light and contaminations prevented during their use.

Do not use reagents over the expiration date.

Signs of reagent deterioration:

- Presence of particles and turbidity.

ADDITIONAL EQUIPMENT

- Spectrophotometer or colorimeter measuring at 582 nm.

- Matched cuvettes 1,0 cm light path.

- General laboratory equipment^(Note 2).

SAMPLES

- Serum or plasma¹: Not haemolyzed. Use only heparin salts as anticoagulants. Stability: 24 hours at 2-8°C or 15 day at -20°C.

PROCEDURE

1. Allow reagents to reach working temperature before using. A proportional variation of the reaction volumes indicated does not change the result.

2. Assay conditions:

Wavelength: 582 nm (570-590)

Cuvette: 1 cm light path

Temperature: 37°C

3. Adjust the instrument to zero with distilled water.

4. Pipette into a cuvette^(Note 3):

	Blank	Standard	Sample
WR (mL)	1,0	1,0	1,0
Distilled water	50	-	-
Standard ^(Note 1) (µL)	--	50	--
Sample (µL)	--	--	50

5. Mix and read the absorbance (A₁) of the sample against the Blank. Add:

R2 (µL)	Blank	Standard	Sample
	50	50	50

6. Mix and incubate for 4-5 min at 37°C.

7. Read the absorbance (A₂) of sample and standard against Blank. The colour is stable for at least 1 hour.

CALCULATIONS

$$\frac{(A_2 - A_1) \text{ Sample} - (A_2 - A_1) \text{ Blank}}{(A_2 - A_1) \text{ Standard} - (A_2 - A_1) \text{ Blank}} \times 100 (\text{Std. conc.}) = \mu\text{g/dL copper in the sample}$$

Conversion factor: µg/dL x 0,1573 = µmol/L.

QUALITY CONTROL

Control sera are recommended to monitor the performance of assay procedures: SPINROL H Normal and Pathologic (Ref. 1002120 and 1002210).

If control values are found outside the defined range, check the instrument, reagents and calibrator for problems.

Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

REFERENCE VALUES

Male 70 - 140 µg/dL ≅ 11,0 - 22,0 µmol/L

Female 80 - 155 µg/dL ≅ 12,6 - 24,4 µmol/L

These values are for orientation purpose; each laboratory should establish its own reference range.

PERFORMANCE CHARACTERISTICS

Measuring range: From *detection limit* of 2,8 µg/dL to *linearity limit* of 848,8 µg/dL.

If the results obtained were greater than linearity limit, dilute the sample 1/10 with NaCl 9 g/L and multiply the result by 10.

Precision:

	Intra-assay (n=20)		Inter-assay (n=20)	
	Mean (µg/dL)	SD	49,7	195,7
SD	1,10	1,43	1,20	2,20
CV (%)	2,23	0,76	2,30	1,20

Sensitivity: 1µg/dL = 0,000767 (A).

Accuracy: Results obtained using SPINREACT reagents (y) did not show systematic differences when compared with other commercial reagents (x).

The results obtained using 80 samples were the following:

Correlation coefficient (r)²: 1,00.

Regression equation: y= 1,03x + 0,25.

The results of the performance characteristics depend on the analyzer used.

INTERFERENCES

No interferences were observed to bilirubin up to 15 mg/dL, hemoglobin up to 0,5 g/dL and triglycerides up to 1000 mg/dL.

A list of drugs and other interfering substances with copper determination has been reported^{3,4}.

NOTES

1. COPPER CAL: Proceed carefully with this product because due its nature it can get contaminated easily.

2. It is recommended to use disposable material. If glassware is used the material should be scrupulously cleaned with hydrochloric acid 1 N and then thoroughly rinsed it with distilled water.

3. Use clean disposable pipette tips for its dispensation.

4. **SPINREACT has instruction sheets for several automatic analyzers. Instructions for many of them are available on request.**

BIBLIOGRAPHY

1. Burtis A et al. Tietz Textbook of Clinical Chemistry, 3rd ed AACC 1999.
2. Kaplan A et al. Clin Chem The C.V. Mosby Co.1984.
3. Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.
4. Young DS. Effects of disease on Clinical Lab. Tests, 4th ed AACC 2001.

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