
**Quantitative determination of urea
IVD**

Store at 2-8°C

PRINCIPLE OF THE METHOD

 Urea in the sample is hydrolyzed enzymatically into ammonia (NH₄⁺) and carbon dioxide (CO₂).

 Ammonia ions formed reacts with α-ketoglutarate in a reaction catalysed by glutamate dehydrogenase (GLDH) with simultaneous oxidation of NADH to NAD⁺:

 The decrease in concentration of NADH, is proportional to urea concentration in the sample¹.

CLINICAL SIGNIFICANCE

Urea is the final result of the metabolism of proteins; it is formed in the liver from their destruction.

 It can appear the urea elevated in blood (uremia) in: diets with excess of proteins, renal diseases, heart failure, gastrointestinal hemorrhage, dehydration or renal obstruction^{1,4,5}.

Clinical diagnosis should not be made on a single test result; it should integrate clinical and other laboratory data.

REAGENTS

R 1 Buffer	TRIS pH 7.8	80 mmol/L
	α-Ketoglutarate	6 mmol/L
	Urease	75000 U/L
R 2 Enzymes	GLDH	60000 U/L
	NADH	0.32 mmol/L

PREPARATION

All the reagents are ready to use.

STORAGE AND STABILITY

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C, protected from light and contaminations prevented during their use.

Do not use reagents over the expiration date.

Signs of reagent deterioration:

- Presence of particles and turbidity.
- Blank absorbance (A) at 340 nm < 1,00.

ADDITIONAL EQUIPMENT

- SPIN640 / SPIN640Plus Autoanalyzer.
- General laboratory equipment^(Note 1).

SAMPLES

 - Serum or heparinized plasma¹: Do not use ammonium salts or fluoride as anticoagulants.

 - Urine¹: Dilute sample 1/50 in distilled water. Mix. Multiply the results by 50 (dilution factor). Preserve urine samples at pH < 4.

Urea is stable at 2-8°C for 5 days.

REFERENCE VALUES^{4,5}

Serum or plasma:

15-45 mg/dL ≅ 2,5-7,5 mmol/L

Urine:

26 – 43 g/24 h ≅ 428-714 mmol/24 h

These values are for orientation purpose; each laboratory should establish its own reference range.

QUALITY CONTROL

Control sera and calibrators are recommended to monitor the performance of assay procedures: SPINCONTROL H Calibrator, SPINCONTROL H Normal and Pathologic (Ref. 1002011, 1002120 and 1002210).

If control values are found outside the defined range, check the instrument, reagents and technique for problems.

Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

BARCODED REAGENTS LOAD MUST BE PRECEDED OF A SPINREACT "DATABASE" COPY INTO THE ANALYZER SOFTWARE. IT IS AVAILABLE UNDER REQUEST TO SPINREACT.

SPIN640 APPLICATION

TEST INFORMATION		REAGENT VOLUME	
Nº	**	Vol. R1	240 µL
Test	UREA	Vol. R2	60 µL
Full Name	Urea	Vol. R3	
Standard nº	1	Vol. R4	
SAMPLE VOLUME		RESULT SETUP	
Vol. Sample Stand.	3 µL	Unit	Decimal 1
Vol. Sample Increas.		mg/dL Inter.	0
Vol. Sample Dec			
REACTION PARAMETERS			
Reac. Type	Fixed Time	Direction	Decrease
Pri. Wave.	340	Reagent Blank	0-0
Sec. Wave.		React. Time	46-53

SPIN640Plus APPLICATION

EDIT PARAMETERS			
Test	UREA	No.	**
Full name	UREA	Print name	UREA
Reac. Type	Fixed Time	Direction	Decrease
Pri. Wave.	340	Sec. Wave.	
Unit	mg/dL	Decimal	0.1
Reagent Blank	0 - 0	React. Time	53 - 60
Vol. Sample	3 µL	R1	240 µL
Increased		R2	60 µL
Decreased		R3	
Sample blank		R4	

The Calibration is stable until **36 days**. After this period the Calibration must be performed again in order to obtain good results.

PERFORMANCE CHARACTERISTICS
Measuring range: From *detection limit* 0,743 mg/dL to *linearity limit* 400 mg/dL.

If the concentration is greater than linearity limit dilute 1/2 the sample with NaCl 9 g/L and multiply the result by 2.

Precision:

	Intra-assay (n=20)		Inter-assay (n=20)	
	Mean (mg/dL)	37,5	120	40,0
SD	1,05	0,92	1,06	2,07
CV (%)	2,79	0,77	2,65	1,65

Sensitivity: 1 mg/dL = 0,00180 A.

Accuracy: Results obtained using SPINREACT reagents (y) did not show systematic differences when compared with other commercial reagent (x).

The results obtained using 50 samples was the following:

 Correlation coefficient (r)²: 0,98209.

Regression equation y= 1,0343x – 1,2105.

The results of the performance characteristics depend on the analyzer used.

NOTES

1. Glassware and distilled water must be free of ammonia and ammonium salts¹.
2. Calibration with the aqueous standard may cause a systematic error in automatic procedures. In these cases, it is recommended to use a serum Calibrator.
3. Use clean disposable pipette tips for its dispensation.

BIBLIOGRAPHY

1. Kaplan A. Urea. Kaplan A et al. Clin Chem The C.V. Mosby Co. St Louis. Toronto. Princeton 1984; 1257-1260 and 437 and 418.
2. Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.
3. Young DS. Effects of disease on Clinical Lab. Tests, 4th ed AACC 2001.
4. Burtis A et al. Tietz Textbook of Clinical Chemistry, 3rd ed AACC 1999.
5. Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed AACC 1995.

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