**Quantitative determination of Lipase IVD**

Store at 2-8°C.

**PRINCIPLE OF THE METHOD**

The pancreatic lipase in presence of colipase, desoxycholate and calcium ions, hydrolyses the substrate 1-2-O-dilauryl-rac-glycero-3-glutaric acid-(6' -methylresorufin)-ester. The sequence of reactions involved in the enzymatic direct lipase determination is the following:

1-2-O-dilauryl-rac-glycero-3-glutaric-(6' -methylresorufin)-ester

$$\xrightarrow{\text{Lipase}} \text{1-2-O-dilauryl-rac-glycerol} + \text{Glutaric-6'-methylresorufin-ester (no}$$

$$\text{stable)} \xrightarrow{\text{OH}^-} \text{Glutaric acid} + \text{Methylresorufin}$$

The rate of methylresorufin formation, measured photometrically, is proportional to the catalytic concentration of lipase present in the sample.

**CLINICAL SIGNIFICANCE**

Lipase (LPS) is a pancreatic enzyme necessary for the absorption and digestion of nutrients that catalyzes the hydrolysis of glycerol esters of fatty acids. Determination of LPS is used for diagnosis of diseases of pancreas such as acute and chronic pancreatitis and obstruction of the pancreatic duct<sup>1,7,8</sup>. Clinical diagnosis should not be made on a single test result; it should integrate clinical and other laboratory data.

**REAGENTS**

<b>R 1</b>	TRIS pH 8.3	40 mmol/L
Buffer	Colipase	≥ 1 mg/L
	Desoxycholate	1,8 mmol/L
	Taurodesoxycholate	≥ 20 mmol/L
<b>R 2</b>	Tartrate pH 4,0	15 mmol/L
Substrate	Lipase Substrate	≥ 0,7 mmol/L
(micro-emulsion)	Calcium chloride (CaCl <sub>2</sub> )	≥ 0,05 mmol/L

**PREPARATION**

The reagents are ready to use.

**STORAGE AND STABILITY**

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C, protected from light and contaminations prevented during their use. Do not use reagents over the expiration date.

**Signs of reagent deterioration:**

- Presence of particles and turbidimetry (note 1)
- Blank absorbance (A) at 570 nm ≥ 1,00.
- R2 is a turbid orange-colored micro-emulsion, discard if turning to red.

**ADDITIONAL EQUIPMENT**

- Spintech 240 Autoanalyzer.
- General laboratory equipment <sup>(note 2)</sup>

**SAMPLES**

Fresh serum samples are recommended.

Spinreact does not recommend the use of EDTA and citrate plasma samples. These sample types could only be used if different reference ranges are established by each laboratory using different sample types other than fresh serum.

Avoid repeated freezing and thawing. Stability: 14 days at 2-8°C and at -20°C.

**QUALITY CONTROL**

Control sera are recommended to monitor the performance of assay procedures: SPINTROL H Normal and Pathologic (Ref. 1002120 and 1002210).

If control values are found outside the defined range, check the instrument, reagents and technique for problems.

Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

**REFERENCE VALUES**

The normal range for pancreatic lipase when measured by the DGGR method is 13-60 U/L<sup>8</sup>.

These values are for orientation purpose; each laboratory should establish its own reference range.

**INTERFERENCES**

Bilirubin up to 75 mg/dL, hemoglobin up to 110 mg/dL and ascorbic acid up to 120 mg/dL do not interfere at the tested Lipase concentrations (40 U/L). 600 mg/dL Intralipid (1300 mg/dL triglycerides) do not interfere.

A list of drugs and other interfering substances with lipase determination has been reported<sup>2, 3</sup>.

**NOTES**

- In some storage conditions (i.e. storage at a temperature lower than the one indicate) a precipitate may appear in the vial, it is not recommended to resuspend the precipitate as it could alter the functionality.
- In order to avoid contamination, it is recommended to use disposable material.

**BARCODED REAGENTS LOAD MUST BE PRECEDED OF A SPINREACT "DATABASE" COPY INTO THE ANALYZER SOFTWARE. IT IS AVAILABLE UNDER REQUEST TO SPINREACT.**

**SPINTECH 240 APPLICATION**

Item Name LIPASE			
<b>DATA INFORMATION</b>		<b>CALIBRATION</b>	
Units	U/L	TYPE	2 Linear Point
Decimals	1		
<b>ANALYSIS</b>		STANDARD	
Type	RATE	#1 0	#4
		#2 CAL	#5
		#3	#6
W.Length 1	570	<b>NORMAL RANGE</b>	
W.Length 2			LOW HIGH
Method	Kinetic	SERUM	MALE FEMALE
<b>CORR</b>		URINE	
SLOPE	INTER		
1.000 x +	0		
Item Name LIPASE			
<b>ASPIRATION</b>		<b>DATA PROCESS</b>	
KIND	Single <input checked="" type="checkbox"/> Double	<b>READ</b>	<b>ABSORBANCE LIMIT</b>
		START END	LOW -3.000
		MAIN 37 45	HIGH 3.000
		SUB	
SAMPLE	VOLUME		
3	μL		
REAGENT 1	300 μL		
REAGENT 2	60 μL		
			<b>ENDPOINT LIMIT</b>
			LINEAR CHECK (%) 90
Third Mix	<input checked="" type="checkbox"/> OFF <input type="checkbox"/> ON	<b>FACTOR</b>	
R1 Blank	<input checked="" type="checkbox"/> Water <input type="checkbox"/> R1-B	Blank Correction	1.000
<b>MONITOR</b>		<b>PROZONE CHECK</b>	
0 LEVEL POINT	1	START END	LIMIT (%)
SPAN	3.000	FIRST	
		SECOND	<input checked="" type="checkbox"/> Low High
		THIRD	<input checked="" type="checkbox"/> Low High

The calibration is stable until **28 days**. After this period the calibration must be performed again in order to obtain good results.

**PERFORMANCE CHARACTERISTICS**

**1. Trueness:** trueness has been performed using Certified Reference Material. Results obtained with Spinreact reagents show an acceptable proportional systematic deviation at low levels compared to the reference standard.

Theoretical conc. (U/L)	Measured conc. (U/L)	Bias (%)
10,04	10,03	-0,03
20,07	20,77	3,49
40,15	42,83	6,68
80,3	83,79	4,34
160,6	163,47	1,77

**2. Precision:** repeatability and reproducibility were evaluated based on the CLSI protocol EP05-A3:

Repeatability (n= 84)		Reproducibility (n= 90)	
Mean (U/L)	CV (%)	Mean (U/L)	CV (%)
46.8	3.3	49.25	5.0
93.46	2.4	103.2	4.3

**3. Analytical sensitivity:** the observed analytical sensitivity at 570nm was 0.378 mAbs/U/L

**4. Limit of detection (LoD):** the Limit of Detection (LoD) for LIPASE-LQ. DGGR is 1.591 U/L, determined based on the CLSI EP17-A2 protocol.

**5. Measuring range:** from *quantification limit* (LoQ) of 4.32 U/L to *linearity limit* of 700 U/L. LoQ determined based on the CLSI EP17-A2 guidelines. And Linearity determined based on CLSI EP06-A 2<sup>nd</sup> Ed.

**6. Method comparison:** results obtained using SPINREACT reagents (y) compared with other commercial reagents (x) did show slight proportional bias (less than 9%). The results obtained using 50 samples were the following:

Correlation coefficient r= 0.994.

Ordinary least-squares fit: y= 1.313 + 1.076x

The results of the correlation may vary between analyzers and on the samples used.

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