

DHG-9053A

Electrical Thermostatic Dry Oven



Scope of application:

It is applicable for drying, baking, dewaxing and sterilization in places like mine enterprise, laboratory, academy, research unit.

Structure:

_The cabinet body is welded by excellent-quality steel plates or stainless steel plates, the work room is made of stainless steel or sprayed by high temperature resistant silver lacquer, the four corners of internal bladder are semi-circular arcs shape, which is easy for cleaning, novel sculpt, good looking and energy saving.

_The temperature controller with computer digital display is accurate and reliable in controlling temperature.

_The hot air cycle system is made of fan with low noise and proper air pipe, with even temperature in the working room.

Specifications:

Model	DHG-9053
Power voltage	~220V 50Hz
Consume power	1000W
Temperature range	+10~300°C
Temperature fluctuation	±1°C
Working size (W*D*H)	415*370*345(mm)
Shelves per chamber	2
Memo	A