

Thalassemia (Desferal) portable infusion pump



Feature

- ❖ Portable, only 200g, compact & light design
- ❖ Injection mode: rate control mode
- ❖ Compatible with all brands of 1ml-50ml syringes
- ❖ AC power supply or 4AA batteries to work, selectable by user
- ❖ Automatic visible, audible and shaking alarms function
- ❖ User-defined bolus rate and volume
- ❖ Voltage 12V, suitable for ambulance
- ❖ With carry bag to take conveniently
- ❖ Memory all parameters automatically

Product Parameter

Applicable syringes	1ml, 2.5ml, 5ml, 10ml, 20ml, 30ml, 50ml disposable syringes 1ml syringe: 0.1- 10.9ml/hr 2.5ml syringe: 0.1-39.6ml/hr 5ml syringe: 0.1- 75.5ml/hr
Flow Rate	10ml syringe: 0.1-99.9ml/hr 20ml syringe: 0.1-99.9ml/hr 30ml syringe: 0.1-99.9ml/hr 50ml syringe: 0.1-99.9ml/hr
Bolus Rate	0.1ml/hr – 99.9ml/hr