

Digital Thermostatic Water Bath



-electric laboratory water bath

-Digital thermostatic water bath is widely used in concentrated solution, distillation, dipping in chemical reagents, drugs and biologics.

-Microcomputer control temperature system, digital display and liner using stainless steel making.

Specifications:

-DIY Supplies: electrical

-Classification: Laboratory Thermostatic Devices

-Model Number: HH-2

-Heating Power: 600W

-Uniformity: <plus or minus 0.5 °C

-Temperature control range: At room temperature to 100 °C plus or minus 0.1 °C

-Working Size: 310*160*120mm

-Power Supply: 220V 50HZ
