

# GEM Premier 3000

## Blood Gas System

### Specifications

Manufacturer:	Instrumentation Laboratory
Series:	GEM Premier
CLIA Complexity:	Moderate
Sample Types:	Whole Blood
Number of Assays:	9
Maximum Throughput:	40 tests/hour
Sample Cycle Time:	85 (sec)
Sample Size Min:	135 $\mu$ L
Sample Input Type:	Cartridge
Storage Capacity Data:	24000
Storage Capacity QC:	24000
Configuration:	Bench Top
Height:	43 cm (17 in)
Width:	31 cm (12 in)
Depth:	31 cm (12 in)
Weight:	13 kg (29 lb)
Power Supply:	110V-220V
Point Of Care:	Yes

