



## ID-Centrifuge 12 S II



---

### Technical characteristics of ID- CENTRIFUGE 12 S II:

Silent centrifuge for 12 ID cards.

Automatic drum balancing control.

Rotation speed 1030 rpm (acceleration 85g), centrifugation time 10 min.

Speed time and other functions are displayed in English, French, German, Italian, Spanish or Portuguese.

Dimensions (width/height/length): 30 cm / 18 cm / 36 cm.

Weight: 6.2 kg.

Power supply requirements: 110 – 240 V / 50 – 60 Hz

---