



SUN-8122 Twelve Channel Interpretive Electrocardiograph with Color Display



Features

- Portable 12 channel ECG machine
- Work modes: Manual/Auto/Arrhythmia analysis modes
- 10 inch colour touch LCD
- Alphabetic keyboard for patient information
- Built in rechargeable Li-ion battery
- Adaptor 110-230V, 50/60Hz AC power supply
- Real-time ECG waveform, 250 patients save/replay
- Built in USB/RS232 interface(ECGNET software Optional)

Technical Specification

- Lead standard 12leads

