

Thermostatic electric water tank DK-420



Thermostatic electric water tank is widely used for distillation, disinfection, concentration and impregnation of chemical or biological products. It is an ideal tool for laboratory personnel in various large and medium sized universities, such as environmental protection, scientific research, health, epidemic prevention petroleum, metallurgy, chemical industry and medical treatment.

Main technical parameters:

Electrical source: 220 +/- 10V 50HZ

Temperature control range: At room temperature to 100 °C plus or minus 0.1 °C

Heating power: 800W

Uniformity: <plus or minus 0.5°C

Work room(mm): 420*180*150
