

TD4A Tabletop Small Low Speed Centrifuge

Widely used in the fields of blood serum, plasma, and urea of hospital clinic and biochemistry.

Features :



- Metal housing is safe and compact. It is of microprocessor control accuracy high.
- Brush-less DC motor, free maintenance, no powder pollution
- Microprocessor control and **LCD or Digital display** Speed/RCF & Time
- There are 10 rates of speed up and down faster and stable.
- Safety door lock alarm inside
- There are five angle rotors are for options
- ISO9001:2008, ISO13485, CE certificates are available as well as calibration report.

Technical Data

Max.Speed	4,000r/min
Max.Capacity	6x50ml
Timer	0~99min
Dimension	485x320x255mm
Speed Accuracy	±20r/min
Max.RCF	2,390xg
Power supply	AC220V, 50/60Hz
Noise	≤45dBA
Net /Gross weight	18/20KG
Temperature range	Non-refrigerated



Matched Rotors

Order No.	Rotor type	Max.Speed (rpm)	Max.Volume (ml)	Max.RCF (xg)
No31605	Angle rotor	4,000	30x7ml	2,250xg
No31601	Angle rotor	4,000	30x5ml	2,250xg
No31604	Angle rotor	4,000	18x10ml	2,250xg
No31607	Angle rotor	4,000	24x10ml	2,200xg
No31606	Angle rotor	4,000	6x50ml	2,100xg
No31602	Angle rotor	4,000	12x15ml	2,150xg